






AIRWAY MANAGEMENT




- 8-10 Oropharyngeal Airways
- 11 Intubation Kit
- 12-13 Nasopharyngeal Airways
- 14-16 Endotracheal Tubes
- 17 Endobronchial Tubes
- 17 Tracheostomy Tubes
- 18-19 ET Connectors
- 20-23 Intubation Aids

Oropharyngeal Airway

BERMAN

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Firm, polyethylene plastic
- Vented septum
- Built-in bite block
- Sets are ideal for crash carts
- Individually wrapped

INDIVIDUAL

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1508-40	Newborn	Pink	4 cm/ 40 mm	10
1-1508-50	Infant	Blue	5 cm/ 50 mm	10
1-1508-60	Small Child	Black	6 cm/ 60 mm	10
1-1508-70	Child	White	7 cm/ 70 mm	10
1-1508-80	Small Adult	Green	8 cm/ 80 mm	10
1-1508-90	Medium Adult	Yellow	9 cm/ 90 mm	10
1-1508-99	Large Adult	Red	10 cm/ 100 mm	10
1-1508-11	X-Large Adult	Orange	11 cm/ 110 mm	10
1-1506-40	Newborn	Pink	4 cm/ 40 mm	50
1-1506-50	Infant	Blue	5 cm/ 50 mm	50
1-1506-60	Small Child	Black	6 cm/ 60 mm	50
1-1506-70	Child	White	7 cm/ 70 mm	50
1-1506-80	Small Adult	Green	8 cm/ 80 mm	50
1-1506-90	Medium Adult	Yellow	9 cm/ 90 mm	50
1-1506-99	Large Adult	Red	10 cm/ 100 mm	50
1-1506-11	X-Large Adult	Orange	11 cm/ 110 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-5008-00	40, 50, 60, 70, 80, 90, 100, 110 mm (1 of each size per set; 10 sets per pack)	10 Set

Oropharyngeal Airway

TRADITIONAL GUEDEL

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Firm, polyethylene plastic
- One piece construction
- Built-in bite block
- Sets are ideal for crash carts
- Bulk packaged

INDIVIDUAL

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1501-60	Small Child	Black	6 cm/ 60 mm	10
1-1501-70	Child	White	7 cm/ 70 mm	10
1-1501-80	Small Adult	Green	8 cm/ 80 mm	10
1-1501-90	Medium Adult	Yellow	9 cm/ 90 mm	10
1-1500-60	Small Child	Black	6 cm/ 60 mm	50
1-1500-70	Child	White	7 cm/ 70 mm	50
1-1500-80	Small Adult	Green	8 cm/ 80 mm	50
1-1500-90	Medium Adult	Yellow	9 cm/ 90 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-1501-00	60, 70, 80, 90 mm (1 of each size per set/pack)	1 Set

Oropharyngeal Airway

COLOR CODED GUEDEL

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Firm, polyethylene plastic
- Built-in bite block
- Sets are ideal for crash carts
- Individually wrapped

INDIVIDUAL

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1504-40	Newborn	Pink	4 cm/ 40 mm	10
1-1504-50	Infant	Blue	5 cm/ 50 mm	10
1-1504-60	Small Child	Black	6 cm/ 60 mm	10
1-1504-70	Child	White	7 cm/ 70 mm	10
1-1504-80	Small Adult	Green	8 cm/ 80 mm	10
1-1504-90	Medium Adult	Yellow	9 cm/ 90 mm	10
1-1504-99	Large Adult	Red	10 cm/ 100 mm	10
1-1504-11	X-Large Adult	Orange	11 cm/ 110 mm	10

1-1505-40	Newborn	Pink	4 cm/ 40 mm	50
1-1505-50	Infant	Blue	5 cm/ 50 mm	50
1-1505-60	Small Child	Black	6 cm/ 60 mm	50
1-1505-70	Child	White	7 cm/ 70 mm	50
1-1505-80	Small Adult	Green	8 cm/ 80 mm	50
1-1505-90	Medium Adult	Yellow	9 cm/ 90 mm	50
1-1505-99	Large Adult	Red	10 cm/ 100 mm	50
1-1505-11	X-Large Adult	Orange	11 cm/ 110 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-1504-00	40, 50, 60, 70, 80, 90, 100, 110 mm (1 of each size per set/pack)	1

Oropharyngeal Airway

PVC GUEDEL

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Soft, clear PVC plastic
- Built-in bite block
- Sets are ideal for crash carts
- Individually wrapped

INDIVIDUAL

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1510-30	Premature	Clear	3 cm/ 30 mm	10
1-1510-20	Newborn	Pink	4 cm/ 40 mm	10
1-1510-10	Infant	Blue	5 cm/ 50 mm	10
1-1510-01	Small Child	Black	6 cm/ 60 mm	10
1-1510-02	Child	White	7 cm/ 70 mm	10
1-1510-03	Small Adult	Green	8 cm/ 80 mm	10
1-1510-04	Medium Adult	Yellow	9 cm/ 90 mm	10
1-1510-05	Large Adult	Red	10 cm/ 100 mm	10
1-1510-06	X-Large Adult	Orange	11 cm/ 110 mm	10

1-1521-30	Premature	Clear	3 cm/ 30 mm	50
1-1521-20	Newborn	Pink	4 cm/ 40 mm	50
1-1521-10	Infant	Blue	5 cm/ 50 mm	50
1-1521-01	Small Child	Black	6 cm/ 60 mm	50
1-1521-02	Child	White	7 cm/ 70 mm	50
1-1521-03	Small Adult	Green	8 cm/ 80 mm	50
1-1521-04	Medium Adult	Yellow	9 cm/ 90 mm	50
1-1521-05	Large Adult	Red	10 cm/ 100 mm	50
1-1521-06	X-Large Adult	Orange	11 cm/ 110 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-1510-98	30, 40, 50, 60, 70, 80, 90, 100, 110 mm (1 of each size per set/pack)	1

Oropharyngeal Airway

WILLIAMS

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Guides fiber optic bronchoscopy placement
- Built-in bite block
- One piece construction

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1503-09	Adult Female	Yellow	9 cm	10
1-1503-10	Adult Male	Pink	10 cm	10

Patient Safety

BITE BLOCK

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox



- Assists with naso-gastric tube placement
- Accommodates up to a 40 Fr tube
- Individually packaged

ITEM	DESCRIPTION	PK
2060	Bite Block	20

Patient Safety

TOOTH PROTECTOR

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Helps prevent broken or chipped teeth during intubation, bronchoscopy, laryngoscopy, tonsillectomy, and tooth extraction
- Minimizes tooth injury

ITEM	DESCRIPTION	PK
9-3041-03	Adult Tooth Protector	50

Patient Safety

JAW SPREADER KIT

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Helps facilitate breathing during convulsions or seizures
- Universal design firmly holds the tongue down while keeping jaw open in both children and adults
- Dependable, durable construction – will not break or splinter even under extreme pressure


ITEM	DESCRIPTION	PK
9-3042-00	Kit contains: Jaw Spreader (fits children and adults); Pediatric Berman Airway (6 cm); Adult Berman Airway (9 cm); Case	1

Intubation

INTUBATION KIT

LATEX FREE

SunMed



Complete intubation kits available with conventional, fiber optic or fiber optic disposable laryngoscope blades

Each kit includes:

- 8 Oropharyngeal color coded Guedel airways (sizes 40 mm-100 mm)
- 1 Penlite chrome plated laryngoscope handle
- 1 Medium chrome plated laryngoscope handle
- 5 Miller laryngoscope blades (sizes 0-4)
- 4 Macintosh laryngoscope blades (sizes 1-4)
- 4 Extra laryngoscope lamps (2 small, 2 large)
- 1 Adult Magill forceps
- 1 Pediatric Magill forceps
- 2 "C" batteries
- 2 "AA" batteries
- 4 Packets of sterile lubricating jelly
- 2 Endotracheal uncuffed tubes (sizes 2.5, 4 mm)
- 8 Endotracheal cuffed tubes (sizes 5-10 mm)
- 4 Nasopharyngeal airways (sizes 20, 24, 28, 32 Fr)
- 1 10 mL syringe
- 1 Endotracheal tube introducer (bougie) with coude tip
- 1 Adult stylette
- 1 Nylon airway management case

ITEM	DESCRIPTION	PK
8-1050-95	Fiber Optic Intubation Kit Includes items above + Reusable FO Laryngoscope Blades	1
8-1050-85	Conventional Intubation Kit Includes items above + Conventional Laryngoscope Blades	1
8-1050-94	Single-Use Fiber Optic Intubation Kit Includes items above + Disposable FO Single Use Laryngoscope Blades	1
8-1050-80	Nylon Case Only	1

Nasopharyngeal Airway

LATEX FREE ROBERTAZZI

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Blunt tip
- Resists collapse or kinking
- Packaged with large size and specification notation for accurate selection
- Sets are ideal for crash carts
- Individually sterile packaged

INDIVIDUAL

CAT#	SIZE	ID	LENGTH	PK
1-5075-20	20 Fr	3.9 mm	114 mm	10
1-5075-22	22 Fr	4.8 mm	114 mm	10
1-5075-24	24 Fr	5.3 mm	120 mm	10
1-5075-26	26 Fr	6.3 mm	121 mm	10
1-5075-28	28 Fr	6.8 mm	132 mm	10
1-5075-30	30 Fr	7.2 mm	139 mm	10
1-5075-32	32 Fr	7.6 mm	146 mm	10
1-5075-34	34 Fr	8.5 mm	157 mm	10
1-5075-36	36 Fr	9.0 mm	170 mm	10

SETS

ITEM	SIZES INCLUDED (1 of each size per set/package)	PK
1-5075-04	20, 24, 28, 32 Fr + sterile lubricating jelly	1 Set
1-5075-05	28, 30, 32, 34, 36 Fr + sterile lubricating jelly	1 Set
1-5075-06	22, 24, 26, 28, 30, 32 Fr + sterile lubricating jelly	1 Set
1-5075-09	20, 22, 24, 26, 28, 30, 32, 34, 36 Fr + sterile lubricating jelly	1 Set

Nasopharyngeal Airway

ADJUSTABLE FLANGE AIRWAY

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Soft, adjustable flange
- Blunt tip
- Resists collapse or kinking
- Individually sterile packaged

INDIVIDUAL

ITEM	SIZE	ID	LENGTH	PK
1-5072-12	12 Fr	3.0 mm	95 mm	10
1-5072-14	14 Fr	3.5 mm	95 mm	10
1-5072-16	16 Fr	4.0 mm	95 mm	10
1-5072-18	18 Fr	4.5 mm	95 mm	10
1-5072-20	20 Fr	5.0 mm	115 mm	10
1-5072-22	22 Fr	5.5 mm	120 mm	10
1-5072-24	24 Fr	6.0 mm	120 mm	10
1-5072-26	26 Fr	6.5 mm	130 mm	10
1-5072-28	28 Fr	7.0 mm	170 mm	10
1-5072-30	30 Fr	7.5 mm	170 mm	10
1-5072-32	32 Fr	8.0 mm	170 mm	10
1-5072-34	34 Fr	8.5 mm	170 mm	10
1-5072-36	36 Fr	9.0 mm	170 mm	10

Nasopharyngeal Airway

LATEX ROBERTAZZI

STERILE | SINGLE USE | DISPOSABLE

SunMed



- Processed to reduce latex protein
- Smooth surface for easy insertion
- Soft, flexible design reduces trauma to nasal passage
- Blunt tip
- Resists collapse or kinking
- Sets are ideal for crash carts
- Individually sterile packaged

INDIVIDUAL

ITEM	SIZE	ID	LENGTH	PK
1-5076-20	20 Fr	3.9 mm	114 mm	10
1-5076-22	22 Fr	5.5 mm	120 mm	10
1-5076-24	24 Fr	6.0 mm	120 mm	10
1-5076-26	26 Fr	6.5 mm	130 mm	10
1-5076-28	28 Fr	7.0 mm	135 mm	10
1-5076-30	30 Fr	7.5 mm	135 mm	10
1-5076-32	32 Fr	8.0 mm	150 mm	10
1-5076-34	34 Fr	8.5 mm	160 mm	10
1-5076-36	36 Fr	9.0 mm	175 mm	10

SETS


ITEM	SIZES INCLUDED (1 of each size per set/package)	PK
1-5076-04	20, 24, 28, 32 Fr + sterile lubricating jelly	1 Set
1-5076-06	22, 24, 26, 28, 30, 32 Fr + sterile lubricating jelly	1 Set
1-5076-09	20, 22, 24, 26, 28, 30, 32, 34, 36 Fr + sterile lubricating jelly	1 Set

Nasopharyngeal Airway

PVC AIRWAY

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Smooth surface for easy insertion
- Beveled tip
- Standard 15 mm Fitting
- Packaged with large size and specification notation for accurate selection
- Sets are ideal for crash carts
- Individually sterile packaged

INDIVIDUAL

ITEM	SIZE	ID	LENGTH	PK
1-5073-12	12 Fr	3.0 mm	75 mm	10
1-5073-14	14 Fr	3.5 mm	90 mm	10
1-5073-16	16 Fr	4.0 mm	105m	10
1-5073-18	18 Fr	4.5 mm	115 mm	10
1-5073-20	20 Fr	5.0 mm	120 mm	10
1-5073-22	22 Fr	5.5 mm	125 mm	10
1-5073-24	24 Fr	6.0 mm	140 mm	10
1-5073-26	26 Fr	6.5 mm	145 mm	10
1-5073-28	28 Fr	7.0 mm	155 mm	10
1-5073-30	30 Fr	7.5 mm	160 mm	10
1-5073-32	32 Fr	8.0 mm	180 mm	10
1-5073-34	34 Fr	8.5 mm	185 mm	10
1-5073-36	36 Fr	9.0 mm	195 mm	10

SETS


ITEM	SIZES INCLUDED (1 of each size per set/package)	PK
1-5073-05	20, 24, 28, 32, 36 Fr + sterile lubricating jelly	1 Set

Endotracheal Tube
MURPHY EYE

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK



- Oral and nasal use
- Clear, flexible PVC
- Smooth Murphy Eye for alternative ventilation port
- Available cuffed, uncuffed, reinforced and with stylette
- Economical
- High volume, low pressure cuff
- Large barrel cuff
- Standard 15 mm ISO connector
- Reinforced version features a curved, stainless steel spiral



CUFF INFLATION BULB

ITEM	DESCRIPTION	PK
4-0101-00	Cuff Inflation Bulb	10




UNCUFFED

ITEM	ID	OD	LENGTH	SIZE	PK
1-7330-20	2.0 mm	2.9 mm	145 mm	8 Fr	10
1-7330-25	2.5 mm	3.5 mm	150 mm	10 Fr	10
1-7330-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7330-35	3.5 mm	4.9 mm	185 mm	14 Fr	10
1-7330-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7330-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7330-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7330-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7330-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7330-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7330-70	7.0 mm	9.6 mm	310 mm	28 Fr	10




CUFFED

ITEM	ID	OD	LENGTH	SIZE	PK
1-7333-25	2.5 mm	3.5 mm	150 mm	10 Fr	10
1-7333-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7333-35	3.5 mm	4.9 mm	188 mm	14 Fr	10
1-7333-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7333-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7333-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7333-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7333-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7333-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7333-70	7.0 mm	9.6 mm	310 mm	28 Fr	10
1-7333-75	7.5 mm	10.2 mm	315 mm	30 Fr	10
1-7333-80	8.0 mm	10.9 mm	330 mm	32 Fr	10
1-7333-85	8.5 mm	11.5 mm	330 mm	34 Fr	10
1-7333-90	9.0 mm	12.1 mm	330 mm	36 Fr	10
1-7333-95	9.5 mm	12.7 mm	330 mm	38 Fr	10
1-7333-10	10.0 mm	13.6 mm	330 mm	40 Fr	10



CUFFED WITH STYLETTE

ITEM	ID	OD	LENGTH	SIZE	PK
1-7343-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7343-35	3.5 mm	4.9 mm	188 mm	14 Fr	10
1-7343-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7343-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7343-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7343-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7343-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7343-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7343-70	7.0 mm	9.6 mm	310 mm	28 Fr	10
1-7343-75	7.5 mm	10.2 mm	315 mm	30 Fr	10
1-7343-80	8.0 mm	10.9 mm	330 mm	32 Fr	10
1-7343-85	8.5 mm	11.5 mm	330 mm	34 Fr	10
1-7343-90	9.0 mm	12.1 mm	330 mm	36 Fr	10
1-7343-95	9.5 mm	12.7 mm	330 mm	38 Fr	10
1-7343-10	10.0 mm	13.6 mm	330 mm	40 Fr	10



REINFORCED CUFFED


ITEM	ID	LENGTH	SIZE	PK
1-7333-35	3.5 mm	195 mm	14 Fr	10
1-7333-40	4.0 mm	215 mm	16 Fr	10
1-7333-45	4.5 mm	240 mm	18 Fr	10
1-7333-50	5.0 mm	256 mm	20 Fr	10
1-7333-55	5.5 mm	286 mm	22 Fr	10
1-7333-60	6.0 mm	295 mm	24 Fr	10
1-7333-65	6.5 mm	301 mm	26 Fr	10
1-7333-70	7.0 mm	315 mm	28 Fr	10
1-7333-75	7.5 mm	325 mm	30 Fr	10
1-7333-80	8.0 mm	337 mm	32 Fr	10
1-7333-85	8.5 mm	337 mm	34 Fr	10
1-7333-90	9.0 mm	337 mm	36 Fr	10
1-7333-95	9.5 mm	337 mm	38 Fr	10
1-7333-10	10.0 mm	337 mm	40 Fr	10

Endotracheal Tube

ORAL PREFORMED (AGT)(RAE)

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Designed to direct tube downward to rest on patient's chin
- Allows circuit to be positioned out of surgical field; ideal for oral and maxillofacial surgery applications
- Polished Murphy Eye
- 15 mm male fitting included
- High volume, low pressure cuff
- Valve with sensitive pilot balloon indicating cuff status
- Radiopaque strip
- Smooth beveled tip for atraumatic introduction

UNCUFFED

ITEM	ID	LENGTH	SIZE	PK
1-7320-30	3.0 mm	160 mm	12 Fr	10
1-7320-35	3.5 mm	180 mm	14 Fr	10
1-7320-40	4.0 mm	200 mm	16 Fr	10
1-7320-45	4.5 mm	220 mm	18 Fr	10
1-7320-50	5.0 mm	240 mm	20 Fr	10
1-7320-55	5.5 mm	270 mm	22 Fr	10

CUFFED

ITEM	ID	LENGTH	SIZE	PK
1-7321-40	4.0 mm	200 mm	16 Fr	10
1-7321-45	4.5 mm	220 mm	18 Fr	10
1-7321-50	5.0 mm	240 mm	20 Fr	10
1-7321-55	5.5 mm	270 mm	22 Fr	10
1-7321-60	6.0 mm	280 mm	24 Fr	10
1-7321-65	6.5 mm	290 mm	26 Fr	10
1-7321-70	7.0 mm	300 mm	28 Fr	10
1-7321-75	7.5 mm	310 mm	30 Fr	10
1-7321-80	8.0 mm	320 mm	32 Fr	10
1-7321-85	8.5 mm	320 mm	34 Fr	10
1-7321-90	9.0 mm	320 mm	36 Fr	10

Endotracheal Tube

NASAL PREFORMED (AGT)(RAE)

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Designed to direct tube over patient's forehead, reducing pressure on nares
- Allows circuit to be positioned out of surgical field; ideal for oral and maxillofacial surgery applications
- Polished Murphy Eye
- 15 mm male fitting included
- High volume, low pressure cuff
- Valve with sensitive pilot balloon indicating cuff status
- Radiopaque strip
- Smooth beveled tip for atraumatic introduction

UNCUFFED

ITEM	ID	LENGTH	SIZE	PK
1-7322-30	3.0 mm	160 mm	12 Fr	10
1-7322-35	3.5 mm	180 mm	14 Fr	10
1-7322-40	4.0 mm	200 mm	16 Fr	10
1-7322-45	4.5 mm	220 mm	18 Fr	10
1-7322-50	5.0 mm	240 mm	20 Fr	10
1-7322-55	5.5 mm	270 mm	22 Fr	10

CUFFED

ITEM	ID	LENGTH	SIZE	PK
1-7323-40	4.0 mm	200 mm	16 Fr	10
1-7323-45	4.5 mm	220 mm	18 Fr	10
1-7323-50	5.0 mm	240 mm	20 Fr	10
1-7323-55	5.5 mm	270 mm	22 Fr	10
1-7323-60	6.0 mm	280 mm	24 Fr	10
1-7323-65	6.5 mm	290 mm	26 Fr	10
1-7323-70	7.0 mm	300 mm	28 Fr	10
1-7323-75	7.5 mm	310 mm	30 Fr	10
1-7323-80	8.0 mm	320 mm	32 Fr	10

Endobronchial Tube

DUAL LUMEN

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Color coded clear (tracheal) and blue (bronchial) for easy differentiation of appropriate lumen, pilot balloon and cuff
- Soft PVC
- High volume, low pressure cuff
- Radiopaque strip
- Graduated depth marks
- Suction catheter included
- Bronchial stylette included
- 15 mm Y-connector included with each tube

LEFT

ITEM	DESCRIPTION	SIZE	PK
1-7370-35	Endobronchial Robertshaw Left	35 Fr	1
1-7370-37	Endobronchial Robertshaw Left	37 Fr	1
1-7370-39	Endobronchial Robertshaw Left	39 Fr	1
1-7370-41	Endobronchial Robertshaw Left	41 Fr	1

RIGHT

ITEM	DESCRIPTION	SIZE	PK
1-7371-35	Endobronchial Robertshaw Right	35 Fr	1
1-7371-37	Endobronchial Robertshaw Right	37 Fr	1
1-7371-39	Endobronchial Robertshaw Right	39 Fr	1
1-7371-41	Endobronchial Robertshaw Right	41 Fr	1

Tracheostomy Tube

TRACHEOSTOMY

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK

SunMed



- Obturator and tube clearly marked
- 15 mm connector swivel adapter
- Soft PVC
- Secure neck band retainer
- Cushion neck band comfortable for patient
- Radiopaque strip embedded for exact location of tube
- Low pressure cuff

UNCUFFED

ITEM	ID	LENGTH	SIZE	PK
1-7390-30	3.0 mm	56.5 mm	12 Fr	10
1-7390-35	3.5 mm	56.5 mm	14 Fr	10
1-7390-40	4.0 mm	56.5 mm	16 Fr	10
1-7390-45	4.5 mm	56.5 mm	18 Fr	10
1-7390-50	5.0 mm	62 mm	20 Fr	10

CUFFED

ITEM	ID	LENGTH	SIZE	PK
1-7393-50	5.0 mm	62 mm	20 Fr	10
1-7393-55	5.5 mm	67 mm	22 Fr	10
1-7393-60	6.0 mm	67 mm	24 Fr	10
1-7393-65	6.5 mm	77 mm	26 Fr	10
1-7393-70	7.0 mm	80 mm	28 Fr	10
1-7393-75	7.5 mm	84 mm	30 Fr	10
1-7393-80	8.0 mm	90 mm	32 Fr	10
1-7393-85	8.5 mm	92 mm	34 Fr	10
1-7393-90	9.0 mm	97 mm	36 Fr	10
1-7393-95	9.5 mm	97 mm	38 Fr	10
1-7393-10	10.0 mm	97 mm	40 Fr	10

ET Tube Connector

RACINE UNIVERSAL CONNECTOR

LATEX FREE | REUSABLE



- Allows movement of breathing circuit at patient end
- Eliminates the need for connectors of varying types and sizes
- Patented diaphragm allows quick and easy tube insertion which provides a firm, tight seal
- Aluminum and silicone construction
- Fits 7 mm to 10.5 mm tubes

ITEM	DESCRIPTION	PK
3-3013-03	4" Racine Universal Connector	1

ET Tube Connector

FLEXIBLE METAL ADAPTER

LATEX FREE | REUSABLE



- 15 mm (M) x 15 mm (F)
- Permits positioning of tracheal tube Y-piece and tube connectors at any convenient angle
- Low 15 mL dead space
- Chrome plated brass

ITEM	DESCRIPTION	PK
3-0186-15	4.75" Flexible Metal Adapter	1

ET Tube Connector

SYNTHETIC RUBBER CONNECTOR

LATEX FREE | REUSABLE



- 15 mm (M) x 15 mm (F)
- Steam autoclavable
- Delrin fitting on both ends
- Lightweight

ITEM	DESCRIPTION	PK
3-0046-03	3" Synthetic Rubber Connector	1
3-0046-06	6" Synthetic Rubber Connector	1
3-0046-09	9" Synthetic Rubber Connector	1
3-0046-13	13" Synthetic Rubber Connector	1



FOR COMPLETE
Customer Service
Technical Information
Care/Cleaning Instructions
Product Literature

CALL
800-433-2797
VISIT
www.Sun-Med.com
EMAIL
info@Sun-Med.com

ET Tube Connector

STAINLESS STEEL CONNECTORS

LATEX FREE | REUSABLE

SunMed



STRAIGHT

ITEM	SIZE	PK
3-3100-40	4 mm x 15 mm (M)	1
3-3100-50	5 mm x 15 mm (M)	1
3-3100-60	6 mm x 15 mm (M)	1
3-3100-70	7 mm x 15 mm (M)	1
3-3100-80	8 mm x 15 mm (M)	1
3-3100-90	9 mm x 15 mm (M)	1



T-ADAPTER (TIB)

- For introducing or sampling gases From the breathing system
- 22 mm (F) x 22 mm (M)
- Male connection features tubing size 0.25" nipple inside

ITEM	DESCRIPTION	PK
3-0196-00	T-Adapter (TIB)	1



ACUTE 60° NO PORT

ITEM	SIZE	PK
3-3101-40	4 mm x 15 mm (M)	1
3-3101-50	5 mm x 15 mm (M)	1
3-3101-60	6 mm x 15 mm (M)	1
3-3101-70	7 mm x 15 mm (M)	1
3-3101-80	8 mm x 15 mm (M)	1
3-3101-90	9 mm x 15 mm (M)	1



ELBOW

- Chrome plated brass Hanks-Rackow elbow
- 15 mm (M) machine end
- Patient end is nominal 22 mm OD x 15 mm ID

ITEM	DESCRIPTION	PK
3-0169-00	Elbow	1



90° NO PORT

ITEM	SIZE	PK
3-3120-40	4 mm x 15 mm (M)	1
3-3120-50	5 mm x 15 mm (M)	1
3-3120-60	6 mm x 15 mm (M)	1
3-3120-70	7 mm x 15 mm (M)	1
3-3120-80	8 mm x 15 mm (M)	1
3-3120-90	9 mm x 15 mm (M)	1



CORRUGATED CONNECTOR

- Stainless steel coupling for 22 mm corrugated tubes
- May be used to connect breathing bags with 22 mm necks to 22 mm corrugated tubes
- 22 mm (M) x 22 mm (M)

ITEM	DESCRIPTION	PK
3-0183-22	Corrugated Connector	1



WASHINGTON T-TUBE


- Similar to Fletcher T-Tube
- 15 mm (F) patient end x 15 mm (M)
- Turns any size corrugated breathing tube (with proper adapters) into an exhalation tube for older children and adults

ITEM	DESCRIPTION	PK
3-0201-15	Washington T-Tube	1

Introducer
BOUGIE

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Pediatric size features calibrated markings on both ends allowing it to be used as a straight or coudé tip
- Pediatric size fits in 4 mm to 6 mm tubes
- Adult size fits in 6 mm to 11 mm tubes

ITEM	DESCRIPTION	SIZE	PK
9-0211-70	Pediatric, Coudé Tip	10 Fr x 70 cm	10
9-0212-70	Adult, Coudé Tip	15 Fr x 70 cm	10
9-0212-72	Adult, Straight Tip	15 Fr x 70 cm	10

Introducer
MALLEABLE INTRODUCER

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Combines benefits of an introducer and stylette into one unit
- Pediatric size features calibrated markings on both ends allowing it to be used as a straight or coudé tip
- Adjustable stop

ITEM	DESCRIPTION	SIZE	PK
9-0213-82	Pediatric Malleable, Coudé & Straight Tip	10 Fr x 70 cm	10
9-0213-92	Adult, Malleable, Coudé Tip	15 Fr x 70 cm	10



FOR COMPLETE
Customer Service
Technical Information
Care/Cleaning Instructions
Product Literature

CALL
800-433-2797
VISIT
www.Sun-Med.com
EMAIL
info@Sun-Med.com



Introducer
BOUGIE-TO-GO

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE



- Innovative compact packaging designed specifically for use by EMS professionals
- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Fits in 6 mm to 11 mm tubes

ITEM	DESCRIPTION	SIZE	PK
9-0212-82	Adult, Coudé Tip	15 Fr x 60 cm	10


Introducer
PORTED INTRODUCER

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE




- Suitable for providing oxygenation during intubation or eliminating light secretions
- Suction and extract tracheal gas samples for analysis
- 70 cm x hose barb fitting included
- 15 mm male connector included
- Marked to determine depth
- Coudé tip with two eyes

ITEM	DESCRIPTION	SIZE	PK
9-0212-80	Adult, Ported	15 Fr x 70 cm	10



Coudé tip with two eyes



Hose barb fitting

Stylette
PVC COATED

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE




- Malleable aluminum
- Sealed in a PVC sheath

ITEM	DESCRIPTION	SIZE	PK
9-0203-06	SunSlip™ with Sheath	10" x 6 Fr	20
9-0203-10	SunSlip™ with Sheath	11.5" x 10 Fr	20
9-0203-14	SunSlip™ with Sheath	12.5" x 14 Fr	20

Stylette
COATED

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE



- Color coded adjustable stopper
- Coated surface reduces the need for lubricants

ITEM	DESCRIPTION	SIZE	PK
9-0209-71	SunGlide™ Stylette	13" x 6 Fr	25
9-0209-72	SunGlide™ Stylette	15.5" x 8 Fr	25
9-0209-73	SunGlide™ Stylette	18" x 10 Fr	25

Stylette
FIBER OPTIC

LATEX FREE | REUSABLE



- Adjustable stop
- Single use stops available

ITEM	DESCRIPTION	SIZE	PK
9-0204-18	Adult Stylette	18"	1
9-0204-14	Pediatric Stylette	14"	1
9-0240-59	Replacement Lamp	--	1
9-0210-05	Replacement Stops	--	5

Stylette
LIGHTED

LATEX FREE | SINGLE USE | DISPOSABLE



- Lighted
- Adjustable stop


ITEM	DESCRIPTION	SIZE	PK
9-0210-01	Adult Stylette	15.75"	1

Stylette

ALUMINUM

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Aluminum
- Malleable

ITEM	DESCRIPTION	SIZE	PK
9-0204-00	Adult Stylette	16" x 10 Fr	1
9-0204-10	Adult Stylette	16" x 10 Fr	10
9-0204-12	Pediatric Stylette	12" x 10 Fr	1
9-0204-22	Pediatric Stylette	12" x 10 Fr	10

Stylette

ALUMINUM

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Malleable with rounded end
- Satin finish to reduce friction
- Pediatric stylette fits in 2-5 mm tubes
- Adult stylette fits in 6-11 mm tubes

ITEM	DESCRIPTION	SIZE	PK
9-0204-11	Pediatric Stylette	11" x 6 Fr	25
9-0204-25	Adult Stylette	16" x 9 Fr	25

Stylette

STAINLESS STEEL

LATEX FREE | REUSABLE

SunMed



- Surgical stainless steel
- Malleable

ITEM	DESCRIPTION	SIZE	PK
9-0207-03	Adult Stylette	16.5" x 6 Fr	1
9-0207-01	Pediatric Stylette	11.5" x 6 Fr	1

Stylette

COPPER

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Solid copper version is malleable with adjustable stop
- Hollow copper tube features suction port

ITEM	DESCRIPTION	SIZE	PK
9-0203-00	Solid Copper Stylette	14.75" x 8 Fr	1
9-0206-01	Hollow Copper Stylette	16.5"	1