



ANESTHESIA



- 46-71 Fiber Optic Laryngoscopes
- 72-92 Conventional Laryngoscopes
- 93-95 Masks
- 96-98 Pressure Infusion
- 99 iGuard
- 100-101 Peripheral Nerve Stimulators
- 102-103 Anesthesia Accessories

GreenLine®/D™ Fiber Optic Laryngoscope
ALL-METAL MACINTOSH

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

SunMed




- Flexible fiber optic bundle protected in black plastic sheath
- Designed with three ball bearings in the heel for secure handle attachment
- Macintosh blades feature beaded tip reduces tissue trauma
- Constructed of surgical grade 303/304 stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5335-00	Neonate	0	80 mm	16 mm	20
5-5335-01	Infant	1	92 mm	22 mm	20
5-5335-02	Child	2	100 mm	21 mm	20
5-5335-03	Medium Adult	3	130 mm	22 mm	20
5-5335-35	Ext. Med. Adult	3.5	145 mm	25 mm	20
5-5335-04	Large Adult	4	155 mm	25 mm	20

GreenLine®/D™ Fiber Optic Laryngoscope
ALL-METAL MILLER

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

SunMed



- Flexible fiber optic bundle protected in black plastic sheath
- Designed with three ball bearings in the heel for secure handle attachment
- Constructed of surgical grade 303/304 stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5336-20	Premature	00	67 mm	10 mm	20
5-5336-00	Neonate	0	77 mm	10 mm	20
5-5336-01	Infant	1	92 mm	10 mm	20
5-5336-15	Small Child	1.5	130 mm	12 mm	20
5-5336-02	Child	2	100 mm	12 mm	20
5-5336-03	Medium Adult	3	130 mm	14 mm	20
5-5336-04	Large Adult	4	155 mm	18 mm	20

GreenLine®/D™ Fiber Optic Laryngoscope ALL-METAL WIS-HIPPLE

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

SunMed



- Flexible fiber optic bundle protected in black plastic sheath
- Designed with three ball bearings in the heel for secure handle attachment
- Constructed of surgical grade 303/304 stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5334-15	Small Child	1.5	115 mm	12 mm	20

GreenLine®/D™ Fiber Optic Laryngoscope ALL-METAL PHILLIPS

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

SunMed



- Flexible fiber optic bundle protected in black plastic sheath
- Designed with three ball bearings in the heel for secure handle attachment
- Macintosh blades feature beaded tip reduces tissue trauma
- Constructed of surgical grade 303/304 stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5339-01	Infant	1	113 mm	11 mm	20
5-5339-02	Child	2	166 mm	17 mm	20

GreenLine®/D™ Fiber Optic Laryngoscope HANDLE WITH GRIP

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

SunMed



5-5338-11 5-5338-10 5-5338-09 5-5338-12


- 5 year sterility and battery life
- Textured grip surface for more secure handling
- One complete, fully-contained unit
- For use with all GreenLine® laryngoscope blades
- LED illumination provides exceptional, clear viewing
- Complies with ISO 7376 Standard
- Mercury-free batteries pre-installed

ITEM	DESCRIPTION	SIZE	PK
5-5338-09	Textured Grip Disp. Handle	Medium	20
5-5338-10	Textured Grip Disp. Handle	Penlite	20
5-5338-11	Textured Grip Disp. Handle	Stubby	20
5-5338-12	Textured Grip Disp. Handle	Standard	20

For mor information, go to:

SunOneSystem.com

BLADE MEASUREMENTS



SunOne™
SYSTEM

PREVENT THE RISK OF CROSS CONTAMINATION with Single Patient Use Fiber Optic Laryngoscope Blades & Handles

GreenLine®/D™ Fiber Optic Laryngoscope
MACINTOSH

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

SunMed



- Disposable blade solves contamination problems and eliminates the cost and time spent cleaning blades and returning them for use
- Answers the clinician's request for a non-plastic, single-use blade suitable for everyday hospital use
- Polished acrylic stem produces exceptional illumination
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5332-01	Infant	1	92 mm	22 mm	20
5-5332-02	Child	2	100 mm	21 mm	20
5-5332-03	Medium Adult	3	130 mm	22 mm	20
5-5332-35	Ext. Med. Adult	3.5	145 mm	22 mm	20
5-5332-04	Large Adult	4	155 mm	25 mm	20

GreenLine®/D™ Fiber Optic Laryngoscope
MILLER


LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

SunMed



- Disposable blade solves contamination problems and eliminates the cost and time spent cleaning blades and returning them for use
- Answers the clinician's request for a non-plastic, single-use blade suitable for everyday hospital use
- Polished acrylic stem produces exceptional illumination
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5333-20	Premature	00	67 mm	10 mm	20
5-5333-00	Newborn	0	77 mm	10 mm	20
5-5333-01	Infant	1	102 mm	10 mm	20
5-5333-02	Child	2	153 mm	12 mm	20
5-5333-03	Medium Adult	3	195 mm	14 mm	20
5-5333-04	Large Adult	4	205 mm	18 mm	20

BLADE MEASUREMENTS


GreenLine®/D™ Fiber Optic Laryngoscope ROBERTSHAW

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

SunMed



- Disposable blade solves contamination problems and eliminate the cost and time spent cleaning blades and returning them for use
- Answers the clinician's request for a non-plastic single-use blade suitable for everyday hospital use
- Polished acrylic stem produces exceptional illumination
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5337-00	Newborn	0	94 mm	11 mm	20
5-5337-01	Infant	1	112 mm	13 mm	20

GreenLine®/D™ Fiber Optic Laryngoscope FLEX TIP

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

SunMed



- Disposable blade solves contamination problems and eliminates the cost and time spent cleaning blades and returning them for use
- Answers the clinician's request for a non-plastic single-use blade suitable for everyday hospital use
- Polished acrylic stem produces exceptional illumination
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5382-03	Medium Adult	3	195 mm	18 mm	10
5-5382-04	Large Adult	4	205 mm	18 mm	10

GreenLine®/D™ Fiber Optic Laryngoscope HANDLE

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

SunMed



- One complete, fully-contained unit
- For use with all GreenLine® laryngoscope blades
- LED illumination provides exceptional, clear viewing
- Ribbed plastic handle for comfortable grip
- Complies with ISO 7376 Standard
- Non-sterile; packaged clean
- Mercury-free batteries pre-installed

ITEM	DESCRIPTION	SIZE	PK
5-0236-69	Ribbed Grip Disp. Handle	Large	20

GreenLine®/D™ Fiber Optic Laryngoscope PENLITE HANDLE

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

SunMed



- Smooth grip for comfortable handling
- One complete, fully-contained unit
- For use with all GreenLine® laryngoscope blades
- LED illumination provides exceptional, clear viewing
- Complies with ISO 7376 Standard
- Non-sterile; packaged clean
- Mercury-free batteries pre-installed

ITEM	DESCRIPTION	SIZE	PK
5-0236-70	Smooth Grip Disp. Handle	Penlite	20

GreenLine® Fiber Optic Laryngoscope
MACINTOSH GERMAN

LATEX FREE | REUSABLE



- Built with an integral, extra large fiber optic bundle delivering a cool, white light
- Design eliminates cracks and crevices allowing superior cleaning and the prevention of cross contamination
- Not chrome plated
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5272-01	Infant	1	93 mm	15 mm	1
5-5272-02	Child	2	100 mm	21 mm	1
5-5272-13	Strong Curve Medium Adult	3	135 mm	20 mm	1
5-5272-03	Medium Adult	3	135 mm	20 mm	1
5-5272-04	Large Adult	4	155 mm	20 mm	1

GreenLine® Fiber Optic Laryngoscope
MACINTOSH GERMAN KIT

LATEX FREE | REUSABLE



- Conveniently packaged laryngoscope kit with handle and blades
- High impact plastic case for easy transport
- Sturdy foam insert carefully protects components
- Includes medium handle with most popular Macintosh German profile blade sizes
- Reusable – all laryngoscope components can be reprocessed

ITEM	DESCRIPTION	PK
5-5272-44	Medium Handle (5-0236-09) Macintosh Blades: 1 (5-5272-01), 2 (5-5272-02), 3 (5-5272-03), 4 (5-5272-04) Replacement Lamp (5-0240--53) Case (8-1050-73)	1

GreenLine® Fiber Optic Laryngoscope
MACINTOSH AMERICAN

LATEX FREE | REUSABLE



- Flange begins at the lock of blade, continues to where the curvature changes direction, and tapers off
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5231-00	Newborn	0	80 mm	15 mm	1
5-5231-01	Infant	1	87 mm	15 mm	1
5-5231-02	Child	2	108 mm	15 mm	1
5-5231-03	Medium Adult	3	128 mm	24 mm	1
5-5231-35	Ext. Med. Adult	3.5	144 mm	22 mm	1
5-5231-00	Large Adult	4	159 mm	15 mm	1

GreenLine® Fiber Optic Laryngoscope
MACINTOSH AMERICAN KIT


LATEX FREE | REUSABLE



- Conveniently packaged laryngoscope kit with handle and blades
- High impact plastic case for easy transport
- Sturdy foam insert carefully protects components
- Includes medium handle with most popular Macintosh American profile blade sizes
- Reusable – all laryngoscope components can be reprocessed


ITEM	DESCRIPTION	PK
5-5231-44	Medium Handle (5-0236-09) Macintosh Blades: 1 (5-5231-01), 2 (5-5231-02), 3 (5-5231-03), 4 (5-5231-04) Replacement Lamp (5-0240--53) Case (8-1050-73)	1

BLADE MEASUREMENTS



GreenLine® Fiber Optic Laryngoscope
IMPROVED VIEW AMERICAN

LATEX FREE | REUSABLE




- Channel extends from heel to mid-blade for improved visualization of the cords
- Channel design helps guide endotracheal tube insertion
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5254-02	Child	2	100 mm	18 mm	1
5-5254-03	Medium Adult	3	130 mm	24 mm	1
5-5254-35	Ext. Med. Adult	3.5	145 mm	25 mm	1
5-5254-04	Large Adult	4	155 mm	26 mm	1

GreenLine® Fiber Optic Laryngoscope
FLANGE-LESS AMERICAN

LATEX FREE | REUSABLE



- Absent upper flange on curved blade is effective for use on patients with limited mouth opening, prominent incisors, receding mandible, short thick neck, or larynx in an extreme anterior anatomical position
- Design greatly reduces the potential trauma during laryngoscopy
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5242-03	Medium Adult	3	130 mm	8 mm	1
5-5242-04	Large Adult	4	158 mm	8 mm	1

GreenLine® Fiber Optic Laryngoscope
BLECHMAN AMERICAN

LATEX FREE | REUSABLE



- Angled tip elevates the epiglottis in supine patient
- Tongue is swept to the left allowing the blade to tip upward for effective viewing from the right side
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5251-03	Medium Adult	3	148 mm	13 mm	1


GreenLine® Fiber Optic Laryngoscope
PRISMVIEW AMERICAN

LATEX FREE | REUSABLE



- Designed with an optically polished prism with 30° refraction without image inversion
- Ideal when the view of the vocal cords is impaired, (especially the most anterior larynx), during nasal intubation, and for post operative examination of the larynx
- Prevent condensation by immersing blade in warm water bath for 10 minutes prior to use
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5229-02	Child	2	100 mm	17 mm	1
5-5229-03	Medium Adult	3	130 mm	24 mm	1
5-5229-35	Ext. Med. Adult	3.5	144 mm	25 mm	1
5-5229-04	Large Adult	4	155 mm	26 mm	1



GreenLine® Fiber Optic Laryngoscope
MACINTOSH ENGLISH

LATEX FREE | REUSABLE



- Designed with a gentle curvature
- Flange begins at the lock of blade and continues to the tip as it narrows
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5232-01	Infant	1	92 mm	22 mm	1
5-5232-02	Child	2	100 mm	21 mm	1
5-5232-03	Medium Adult	3	130 mm	22 mm	1
5-5232-04	Large Adult	4	155 mm	25 mm	1
5-5232-05	Extra Lg. Adult	5	175 mm	22 mm	1
5-5243-02	Left Hand Adult	3	130 mm	22 mm	1

GreenLine® Fiber Optic Laryngoscope
POLIO ENGLISH

LATEX FREE | REUSABLE



- Modified Macintosh English profile blade for indirect laryngoscopy
- Blade is offset from handle at an obtuse angle to allow intubation of patients on respirators, in body jackets, as well as other difficult intubation scenarios
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5247-03	Medium Adult Angled at 120°	3	130 mm	23 mm	1

GreenLine® Fiber Optic Laryngoscope

IMPROVED VISION ENGLISH

LATEX FREE | REUSABLE

SunMed



- Channel extends from heel to mid-blade for improved visualization of the cords
- Channel design helps guide endotracheal tube insertion
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

GreenLine® Fiber Optic Laryngoscope

PRISMVIEW ENGLISH

LATEX FREE | REUSABLE

SunMed




- Designed with an optically polished prism with 30° refraction without image inversion
- Ideal when the view of the vocal cords is impaired, (especially the most anterior larynx), during nasal intubation, and for post operative examination of the larynx
- Prevent condensation by immersing blade in warm water bath for 10 minutes prior to use
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5246-02	Child	2	100 mm	21 mm	1
5-5246-03	Medium Adult	3	130 mm	23 mm	1
5-5246-35	Ext. Med. Adult	3.5	145 mm	24 mm	1
5-5246-04	Large Adult	4	155 mm	25 mm	1

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5228-02	Child	2	100 mm	17 mm	1
5-5228-03	Medium Adult	3	130 mm	24 mm	1
5-5228-35	Ext. Med. Adult	3.5	144 mm	25 mm	1
5-5228-04	Large Adult	4	155 mm	26 mm	1

BLADE MEASUREMENTS



GreenLine® Fiber Optic Laryngoscope
SUN-FLEX™ GERMAN

LATEX FREE | REUSABLE

SunMed




- Built with an integral, extra large fiber optic bundle delivering a cool, white light
- Design eliminates cracks and crevices, allowing for superior cleaning and the prevention of cross contamination
- Designed to enable precise control to elevate the epiglottis
- Non-pinching tip
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5282-02	Child	2	102 mm	20 mm	1
5-5282-03	Medium Adult	3	132 mm	23 mm	1
5-5282-04	Large Adult	4	155 mm	26 mm	1

GreenLine® Fiber Optic Laryngoscope
REDUCED FLANGE E-MAC™

LATEX FREE | REUSABLE


SunMed



- Channel extends from heel to mid-blade for improved visualization of the cords
- Channel design helps guide endotracheal tube insertion
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5271-02	Child	2	100 mm	17 mm	1
5-5271-03	Medium Adult	3	130 mm	24 mm	1
5-5271-35	Ext. Med. Adult	3.5	145 mm	24 mm	1
5-5271-04	Large Adult	4	155 mm	26 mm	1

BLADE MEASUREMENTS



GreenLine® Fiber Optic Laryngoscope SIKER MIRROR ENGLISH

LATEX FREE | REUSABLE

SunMed




- Ideal for use on patients with anatomical variations which make intubation with conventional blades difficult or impossible
- The basic portions of the blade make an angle of 135° with a mirror located at flange angle facing spatula for better visualization during intubation
- Distal portion is 3"
- Since mirror inverts the reflected image, operator should exercise caution until familiar with its use
- Curved stylette is recommended for use with the Siker blade
- Immersing blade in warm water bath for 10 minutes prior to use will prevent mirror fogging
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5245-04	Large Adult	4	158 mm	18 mm	1

GreenLine® Fiber Optic Laryngoscope MILLER GERMAN

LATEX FREE | REUSABLE

SunMed



- Built with an integral, extra large fiber optic bundle delivering a cool, white light
- Design eliminates cracks and crevices allowing superior cleaning and the prevention of cross contamination
- Not chrome plated
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5273-20	Premature	00	65 mm	11 mm	1
5-5273-00	Newborn	0	77 mm	12 mm	1
5-5273-01	Infant	1	110 mm	12 mm	1
5-5273-15	Small Child	1.5	130 mm	14 mm	1
5-5273-02	Child	2	152 mm	14 mm	1
5-5273-03	Medium Adult	3	192 mm	15 mm	1
5-5273-04	Large Adult	4	207 mm	16 mm	1
5-5243-02	Left Hand Adult	3	130 mm	22 mm	1

GreenLine® Fiber Optic Laryngoscope
MILLER AMERICAN

LATEX FREE | REUSABLE

SunMed



- Light-on-right American profile
- Miller design is predominant choice among straight blades
- Reduced flange helps minimize trauma
- Distal end curve improves lifting of epiglottis
- Designed to facilitate greater exposure of larynx
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5230-02	Child	2	155 mm	13 mm	1
5-5230-03	Medium Adult	3	195 mm	13 mm	1

GreenLine® Fiber Optic Laryngoscope
SNOW AMERICAN

LATEX FREE | REUSABLE

SunMed




- Combined Miller spatula and the Wis-Foregger flange design can reduce trauma during intubation, especially for use with a cuffed tube
- Curves back from the distal end for a convenient visual field
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5249-03	Medium Adult	3	162 mm	16 mm	1

GreenLine® Fiber Optic Laryngoscope CRANWALL AMERICAN

LATEX FREE | REUSABLE

SunMed



- Reduced flange allows easier insertion through a restricted opening, such as arthritic jaw or facial burns
- Design helps decrease potential damage to upper teeth
- Flange removed directly below the light carrier
- Can easily be used as a left-handed straight blade
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

GreenLine® Fiber Optic Laryngoscope MILLER-PORT AMERICAN

LATEX FREE | REUSABLE

SunMed




- Specifically designed to perform laryngoscopy and intubation with reduced risk of hypoxia
- Integrated tube permits delivery of oxygen or other gas mixtures during intubation
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5236-00	Newborn	0	75 mm	10 mm	1
5-5236-01	Infant	1	100 mm	11 mm	1
5-5236-02	Child	2	157 mm	11 mm	1
5-5236-03	Medium Adult	3	195 mm	12 mm	1


ITEM	DESCRIPTION	SIZE	A	B	PK
5-5248-00	Newborn	0	75 mm	10 mm	1
5-5248-01	Infant	1	102 mm	10 mm	1

BLADE MEASUREMENTS



GreenLine® Fiber Optic Laryngoscope
WISCONSIN AMERICAN

LATEX FREE | REUSABLE




- Straight spatula design with flange that expands slightly toward the distal portion of the blade increases visual field helping to prevent potential trauma during intubation
- Distal portion of the blade is wider and formed slightly to the right, making it well-adapted to lifting the epiglottis
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5234-02	Child	2	135 mm	18 mm	1
5-5234-03	Medium Adult	3	162 mm	17 mm	1
5-5234-04	Large Adult	4	199 mm	21 mm	1

GreenLine® Fiber Optic Laryngoscope
WIS-HIPPLE AMERICAN

LATEX FREE | REUSABLE



- Wisconsin-modified blade with a straight spatula and a large, circular flange
- Flange tends to be straighter than a Wisconsin blade and runs parallel with the spatula
- Designed primarily for use in infants
- Wider, flat tip is well-adapted to lifting the epiglottis
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5234-20	Premature	00	65 mm	13 mm	1
5-5234-00	Newborn	0	76 mm	13 mm	1
5-5234-01	Infant	1	102 mm	13 mm	1
5-5234-15	Small Child	1.5	115 mm	14 mm	1

GreenLine® Fiber Optic Laryngoscope
WHITEHEAD AMERICAN

LATEX FREE | REUSABLE




- Reduced flange increases visualization to facilitate intubation
- Designed to reduce potential damage to incisors
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5244-03	Medium Adult	3	162 mm	15 mm	1

GreenLine® Fiber Optic Laryngoscope
PHILLIPS AMERICAN


LATEX FREE | REUSABLE



- Combines the straight Jackson blade design with the curved distal tip Miller design to allow greater visibility and a 'direct line' approach to the trachea
- Low profile flange reduces the risks of oral damage
- Unique lamp mounting enables deep illumination both downward and inward
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5235-01	Child	1	113 mm	13 mm	1
5-5235-02	Adult	1	166 mm	18 mm	1

BLADE MEASUREMENTS




GreenLine® Fiber Optic Laryngoscope

MILLER ENGLISH

LATEX FREE | REUSABLE

SunMed



- Miller design is the predominant choice among straight blades
- Reduced flange helps minimize trauma
- Distal end curve improves lifting of epiglottis
- Designed to facilitate greater exposure of larynx
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5233-20	Premature	00	67 mm	10 mm	1
5-5233-00	Newborn	0	77 mm	10 mm	1
5-5233-01	Infant	1	102 mm	10 mm	1
5-5233-15	Small Child	1.5	130 mm	12 mm	1
5-5233-02	Child	2	153 mm	12 mm	1
5-5233-03	Medium Adult	3	195 mm	12 mm	1
5-5233-04	Large Adult	4	205 mm	16 mm	1

GreenLine® Fiber Optic Laryngoscope

SEWARD ENGLISH

LATEX FREE | REUSABLE


SunMed



- Originally produced as an infant blade for use in asphyxia neonatorum, it is now available to fit older children
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5239-01	Infant	1	105 mm	9 mm	1
5-5239-02	Child	2	130 mm	16 mm	1


BLADE MEASUREMENTS



GreenLine® Fiber Optic Laryngoscope ROBERTSHAW ENGLISH

LATEX FREE | REUSABLE

SunMed



- Designed specifically for infants and children
- Spatula is gently curved over the distal third and can indirectly lift the epiglottis similar to a Macintosh blade
- Permits binocular vision allowing better judgment of depth and consequently less risk of trauma
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5237-00	Newborn	0	94 mm	11 mm	1
5-5237-01	Infant	1	112 mm	13 mm	1

GreenLine® Fiber Optic Laryngoscope SOPER ENGLISH

LATEX FREE | REUSABLE

SunMed



- Broad flat blade makes it easier to restrict premature and newborn infant tongue movement
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	PK	A	B
5-5240-01	Infant	1	1	100 mm	17 mm

GreenLine® Fiber Optic Laryngoscope OXFORD ENGLISH

LATEX FREE | REUSABLE

SunMed



- Suitable for premature infants, babies and young children up to 4yrs
- Overhang on the open side prevents lips from obscuring vision
- Broad, flat surface aids pediatric patients with extreme cleft palates
- Soft, matte finish virtually eliminates reflection and glare
- Compliant with ISO standard 7376
- Green System compatible
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5238-01	Infant	1	102 mm	16 mm	1

GreenLine® Fiber Optic Laryngoscope

REPLACEMENT HANDLE LAMPS

LATEX FREE | REUSABLE

SunMed



5-0240-52



5-0240-53



5-0240-59



5-0238-60

- For use in fiber optic reusable handles
- Meets stringent manufacturing controls
- Contact Customer Care for correct lamp selection

NEWLY DESIGNED REMOVABLE LAMP CARTRIDGES:




5-0240-61



5-0240-62



5-0240-63



Bulb
Lamp Cartridge
Cartridge Base

ITEM	LAMP TYPE	DESCRIPTION	COMPATIBILITY	PK
5-0240-61	Halogen Xenon	Blue Cartridge Base	GreenLine FO Handles (2015 and newer): 5-0236-09, 5-0236-10, 5-0236-11, 5-0236-15, 5-0236-30, 5-0236-31, 5-0236-35, 5-0236-36, 5-0236-38, 5-0236-39	1
5-0240-62	Eco LED	Orange Cartridge Base	GreenLine SunBrite FO Handles (2015 and newer): 5-0236-90, 5-0236-91, 5-0236-99	1
5-0240-63	High End LED	Yellow Cartridge Base	LED FO Handles (2015 and newer): 5-0236-80, 5-0236-89	1
5-0240-52		Red Cartridge Base	GreenLine FO Handles (pre-2012): 5-0236-09, 5-0236-10, 5-0236-11, 5-0236-15, 5-0236-30, 5-0236-31, 5-0236-35, 5-0236-36, 5-0236-38, 5-0236-39	1
5-0240-53		Gold Cartridge	GreenLine FO Handles (2012-2015): 5-0236-09, 5-0236-10, 5-0236-11, 5-0236-15, 5-0236-30, 5-0236-31, 5-0236-35, 5-0236-36, 5-0236-38, 5-0236-39	1
5-0238-59		No Threads		1
5-0238-60		Large Thread		1
5-0240-59		Black Cartridge Base	GreenLine FO Handles (pre-2007): 5-0236-09, 5-0236-10, 5-0236-11, 5-0236-15, 5-0236-30, 5-0236-31, 5-0236-35, 5-0236-36, 5-0236-38, 5-0236-39	1

GreenLine® Fiber Optic Laryngoscope

CHROME PLATED

LATEX FREE | REUSABLE



- Compatible with all Green System fiber optic blades
- Chrome plated brass
- Bright Halogenated Xenon lamp
- Batteries not included

ITEM	DESCRIPTION	SIZE	PK	BATTERIES
5-0236-09	Chrome Plated Handle	Medium	1	2 "C"
5-0236-10	Chrome Plated Handle	Penlite	1	2 "AA"
5-0236-11	Chrome Plated Handle	Stubby	1	2 "AA"
5-0236-15	Chrome Plated Handle	Large	1	2 "D"

GreenLine® Fiber Optic Laryngoscope

PREMIUM CHROME PLATED

LATEX FREE | REUSABLE



- Compatible with all Green System fiber optic blades
- Chrome plated brass
- Premium LED handles have a lifetime warranty
- Low power use prolongs battery life
- Compliant with ISO standard 7376
- Autoclavable
- Batteries not included

ITEM	DESCRIPTION	SIZE	PK	BATTERIES
5-0236-80	Lifetime LED	Penlite	1	2 "AA"
5-0236-89	Lifetime LED	Medium	1	2 "C"

GreenLine® Fiber Optic Laryngoscope
STAINLESS STEEL

LATEX FREE | REUSABLE



- Compatible with all Green System fiber optic blades
- Bright Halogenated Xenon lamp
- Batteries not included

ITEM	DESCRIPTION	SIZE	PK	BATTERIES
5-0236-30	Stainless Steel	Penlite	1	2 "AA"
5-0236-39	Stainless Steel	Medium	1	2 "C"
5-0236-31	Stainless Steel	Stubby	1	2 "AA"
5-0236-35	Stainless Steel	Large	1	2 "D"
5-0236-36	Stainless Steel	Micro	1	Lithium 123
5-0236-38	Stainless Steel	Mini	1	Lithium 123

GreenLine® Fiber Optic Laryngoscope
PREMIUM STAINLESS STEEL

LATEX FREE | REUSABLE



- Compatible with all Green System fiber optic blades
- Premium LED handles have a lifetime warranty
- Low power use prolongs battery life
- Compliant with ISO standard 7376
- Autoclavable
- Batteries not included

ITEM	DESCRIPTION	SIZE	PK	BATTERIES
5-0246-80	Premium Stainless	Penlite	1	2 "AA"
5-0246-81	Premium Stainless	Stubby	1	2 "AA"
5-0246-86	Premium Stainless	Micro	1	Lithium 123
5-0246-88	Premium Stainless	Mini	1	Lithium 123
5-0246-89	Premium Stainless	Medium	1	2 "C"

GreenLine® Fiber Optic Laryngoscope
SUNBRITE LED

LATEX FREE | REUSABLE



- Compatible with all Green System fiber optic blades
- Chrome plated brass
- Low power use prolongs battery life
- Cool, white illumination
- Autoclavable
- Batteries not included

ITEM	DESCRIPTION	SIZE	PK	BATTERIES
5-0236-90	SunBrite LED	Penlite	1	2 "AA"
5-0236-91	SunBrite LED	Stubby	1	2 "AA"
5-0236-99	SunBrite LED	Medium	1	2 "C"

GreenLine® Fiber Optic Laryngoscope
HOWLAND LOCK

LATEX FREE | REUSABLE



- Natural lifting simplifies laryngeal exposure, making intubation possible even in the most difficult cases
- Built-in leverage decreases prying and reduces possibility of broken teeth
- Compatible with Green System fiber optic handles and blades
- Can often make an otherwise certain tracheotomy unnecessary

ITEM	DESCRIPTION	PK
5-0236-16	Fiber Optic Howland Lock	1

Intubation
INTUBATION KIT

LATEX FREE

SunMed



Not all items in kit shown.

Complete intubation kits available with conventional, fiber optic or fiber optic disposable laryngoscope blades

Each kit includes:

- 8 Oropharyngeal color coded Guedel airways (sizes 40 mm-100 mm)
 - 1 Penlite chrome plated laryngoscope handle
 - 1 Medium chrome plated laryngoscope handle
 - 5 Miller laryngoscope blades (sizes 0-4)
 - 4 Macintosh laryngoscope blades (sizes 1-4)
 - 4 Extra laryngoscope lamps (2 small, 2 large)
 - 1 Adult Magill forceps
 - 1 Pediatric Magill forceps
 - 2 "C" batteries
 - 2 "AA" batteries
 - 4 Packets of sterile lubricating jelly
 - 2 Endotracheal uncuffed tubes (sizes 2.5, 4 mm)
 - 8 Endotracheal cuffed tubes (sizes 5-10 mm)
- 4 Nasopharyngeal airways (sizes 20, 24, 28, 32 Fr)
 - 1 10 mL syringe
 - 1 Endotracheal tube introducer (bougie) with coudé tip
 - 1 Adult stylette
 - 1 Nylon airway management case

ITEM	DESCRIPTION	PK
8-1050-95	Fiber Optic Intubation Kit Includes items above + Reusable FO Laryngoscope Blades	1
8-1050-85	Conventional Intubation Kit Includes items above + Conventional Laryngoscope Blades	1
8-1050-94	Single-Use Fiber Optic Intubation Kit Includes items above + Disposable FO Single Use Laryngoscope Blades	1
8-1050-80	Nylon Case Only	1

GreenLine® Fiber Optic Laryngoscope

DISPOSABLE LARYNGOSCOPE KITS







LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Ideal for ambulances, paramedics and first responders
- Variety of blade/handle combinations in most popular sizes
- All blades are single patient use and disposable
- Available with either single use or reusable handles
- Components arrive in original packaging inside sturdy canvas pouch

ITEM	KIT INCLUDES	PK
5-5332-55	2 Disposable Standard Handles (5-5338-12) Macintosh Blades: 2 (5-5332-02), 3 (5-5332-03), 4 (5-5332-04) Miller Blades: 2 (5-5333-02), 3 (5-5333-03) Canvas Pouch	1
5-5332-57	Disposable Standard Handle (5-5338-12) Disposable Stubby Handle (5-5388-11) Macintosh Blades: 2 (5-5332-02), 3 (5-5332-03), 4 (5-5332-04) Miller Blades: 2 (5-5333-02), 3 (5-5333-03) Canvas Pouch	1
5-5336-54	2 Disposable Penlite Handles (5-5338-10) All-Metal Macintosh Blade: 1 (5-5335-01) All-Metal Miller Blades: 00 (5-5336-20), 0 (5-5336-00), 1 (5-5336-01) Canvas Pouch	1
5-5333-58	2 Disposable Standard Handles (5-5338-12) Macintosh Blades: 2 (5-5332-02), 3 (5-5332-03), 4 (5-5332-04) Miller Blades: 0 (5-5333-20), 1 (5-5333-01), 2 (5-5333-02), 3 (5-5333-03) Canvas Pouch	1
5-5333-57	Chrome Plated Med. Handle (5-0236-09) Macintosh Blades: 2 (5-5332-02), 3 (5-5332-03), 4 (5-5332-04) Miller Blades: 0 (5-5333-20), 1 (5-5333-01), 2 (5-5333-02), 3 (5-5333-03) Canvas Pouch	1
5-5335-47	Chrome Plated Med. Handle (5-0236-09) All-Metal Macintosh Blade: 2 (5-5335-02), 3 (5-5335-03), 4 (5-5335-04) All-Metal Miller Blades: 0 (5-5336-00), 1 (5-5336-01), 2 (5-5336-02), 3 (5336-03) Canvas Pouch	1



GreenLine® Fiber Optic Laryngoscope

REUSABLE LARYNGOSCOPE KITS

LATEX FREE | REUSABLE

SunMed

- All blades and handles compatible with Welch Allyn, Propper, Rusch, Heine and AMS
- Blades available in Macintosh and Miller, American, English and German profiles
- High impact plastic case for easy transport



Macintosh English Profile



Miller English Profile

MACINTOSH ENGLISH PROFILE

ITEM	KIT INCLUDES	PK
5-5232-44	Medium Handle (5-0236-09) Replacement Lamp (5-0240-53) Macintosh Blades: 1 (5-5232-01), 2 (5-5232-02), 3 (5-5232-03), 4 (5-5232-04) Case (8-1050-73)	1 Kit
5-5232-45	Medium Premium LED Handle (5-0236-89) Macintosh Blades: 1 (5-5232-01), 2 (5-5232-02), 3 (5-5232-03), 4 (5-5232-04) Case (8-1050-73)	1 Kit
5-5232-46	Medium SunBrite Handle (5-0236-99) Macintosh Blades: 1 (5-5232-01), 2 (5-5232-02), 3 (5-5232-03), 4 (5-5232-04) Case (8-1050-73)	1 Kit

MILLER ENGLISH PROFILE

ITEM	KIT INCLUDES	PK
5-5233-44	Medium SunBrite Handle (5-0236-99) Miller Blades: 0 (5-5233-00), 1 (5-5233-01), 2 (5-5233-02), 3 (5-5233-03), 4 (5-5233-04) Case (8-1050-71)	1 Kit
5-5233-45	Medium Handle (5-0236-09) Replacement Lamp (5-0240-53) Miller Blades: 0 (5-5233-00), 1 (5-5233-01), 2 (5-5233-02), 3 (5-5233-03), 4 (5-5233-04) Case (8-1050-71)	1 Kit
5-5232-46	Medium Premium LED Handle (5-0236-89) Miller Blades: 0 (5-5233-00), 1 (5-5233-01), 2 (5-5233-02), 3 (5-5233-03), 4 (5-5233-04) Case (8-1050-71)	1 Kit

MACINTOSH AMERICAN PROFILE

ITEM	KIT INCLUDES	PK
5-5231-44	Medium Handle (5-0236-09) Replacement Lamp (5-0240-53) Macintosh Blades: 1 (5-5231-01), 2 (5-5231-02), 3 (5-5231-03), 4 (5-5231-04) Case (8-1050-73)	1 Kit
5-5231-45	Medium Premium LED Handle (5-0236-89) Macintosh Blades: 1 (5-5231-01), 2 (5-5231-02), 3 (5-5231-03), 4 (5-5231-04) Case (8-1050-73)	1 Kit
5-5231-46	Medium SunBrite Handle (5-0236-99) Macintosh Blades: 1 (5-5231-01), 2 (5-5231-02), 3 (5-5231-03), 4 (5-5231-04) Case (8-1050-73)	1 Kit

MACINTOSH / MILLER

ITEM	KIT INCLUDES	PK
5-5233-47	Medium Handle (5-0236-09) Replacement Lamp (5-0240-53) Macintosh Blades: 2 (5-5232-02), 3 (5-5232-03), 4 (5-5232-04) Miller Blades: 0 (5-5233-00), 1 (5-5233-01), 2 (5-5233-02), 3 (5-5233-03) Case (8-1050-77)	1 Kit
5-5233-48	Medium Premium LED Handle (5-0236-89) Macintosh Blades: 2 (5-5232-02), 3 (5-5232-03), 4 (5-5232-04) Miller Blades: 0 (5-5233-00), 1 (5-5233-01), 2 (5-5233-02), 3 (5-5233-03) Case (8-1050-77)	1 Kit
5-5233-49	Medium SunBrite Handle (5-0236-99) Macintosh Blades: 2 (5-5232-02), 3 (5-5232-03), 4 (5-5232-04) Miller Blades: 0 (5-5233-00), 1 (5-5233-01), 2 (5-5233-02), 3 (5-5233-03) Case (8-1050-77)	1 Kit



Macintosh/Miller



Macintosh German Profile



Macintosh/Miller German Profile



SunFlex

MACINTOSH GERMAN PROFILE

ITEM	KIT INCLUDES	PK
5-5272-44	Medium Handle (5-0236-09) Replacement Lamp (5-0240-53) Macintosh Blades: 1 (5-5272-01), 2 (5-5272-02), 3 (5-5272-03), 4 (5-5272-04) Case (8-1050-73)	1 Kit
5-5272-45	Medium Premium LED Handle (5-0236-89) Macintosh Blades: 1 (5-5272-01), 2 (5-5272-02), 3 (5-5272-03), 4 (5-5272-04) Case (8-1050-73)	1 Kit
5-5272-46	Medium SunBrite Handle (5-0236-99) Macintosh Blades: 1 (5-5272-01), 2 (5-5272-02), 3 (5-5272-03), 4 (5-5272-04) Case (8-1050-73)	1 Kit

MACINTOSH / MILLER GERMAN PROFILE

ITEM	KIT INCLUDES	PK
5-5272-47	Medium Handle (5-0236-09) Replacement Lamp (5-0240-53) Macintosh Blades: 2 (5-5272-02), 3 (5-5272-03), 4 (5-5272-04) Miller Blades: 0 (5-5273-00), 1 (5-5273-01), 2 (5-5273-02), 3 (5-5273-03) Case (8-1050-77)	1 Kit
5-5272-48	Medium Premium LED Handle (5-0236-89) Macintosh Blades: 2 (5-5272-02), 3 (5-5272-03), 4 (5-5272-04) Miller Blades: 0 (5-5273-00), 1 (5-5273-01), 2 (5-5273-02), 3 (5-5273-03) Case (8-1050-77)	1 Kit
5-5272-49	Medium SunBrite Handle (5-0236-99) Macintosh Blades: 2 (5-5272-02), 3 (5-5272-03), 4 (5-5272-04) Miller Blades: 0 (5-5273-00), 1 (5-5273-01), 2 (5-5273-02), 3 (5-5273-03) Case (8-1050-77)	1 Kit

MILLER GERMAN PROFILE

ITEM	KIT INCLUDES	PK
5-5273-45	Medium Handle (5-0236-09) Replacement Lamp (5-0240-53) Miller Blades: 0 (5-5273-00), 1 (5-5273-01), 2 (5-5273-02), 3 (5-5273-03), 4 (5-5273-04) Case (8-1050-71)	1 Kit
5-5273-46	Medium Premium LED Handle (5-0236-89) Miller Blades: 0 (5-5273-00), 1 (5-5273-01), 2 (5-5273-02), 3 (5-5273-03), 4 (5-5273-04) Case (8-1050-71)	1 Kit
5-5273-47	Medium SunBrite Handle (5-0236-99) Miller Blades: 0 (5-5273-00), 1 (5-5273-01), 2 (5-5273-02), 3 (5-5273-03), 4 (5-5273-04) Case (8-1050-71)	1 Kit

SUN-FLEX MACINTOSH ENGLISH

ITEM	KIT INCLUDES	PK
5-5282-42	Medium Handle (5-0236-09) Replacement Lamp (5-0240-53) SunFlex Blades: 3 (5282-03), 4 (5-5282-04) Case (8-1050-72)	1 Kit
5-5282-43	Medium Premium LED Handle (5-0236-89) SunFlex Blades: 3 (5282-03), 4 (5-5282-04) Case (8-1050-72)	1 Kit
5-5282-44	Medium SunBrite Handle (5-0236-99) SunFlex Blades: 3 (5282-03), 4 (5-5282-04) Case (8-1050-72)	1 Kit

Conventional Laryngoscope

DISPOSABLE LED MACINTOSH

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

SunMed



- Non-plastic, stainless steel single-use blade
- Eliminates the cost associated with cleaning reusable blades
- Helps prevent cross-contamination
- Patented safety heel inhibits blade tip from contacting handle
- Clean cool white LED illumination
- Low power use ensures long battery life
- Convenient for hospital and EMS use
- Compatible with all conventional ISO system handles
- Individually packaged; non-sterile

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5132-01	Infant	1	90 mm	18 mm	20
5-5132-02	Child	2	100 mm	20 mm	20
5-5132-03	Medium Adult	3	135 mm	21 mm	20
5-5132-04	Large Adult	4	155 mm	23 mm	20

Conventional Laryngoscope

LED STANDARD/D MILLER

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE


SunMed



- Non-plastic, stainless steel single-use blade
- Eliminates the cost associated with cleaning reusable blades
- Helps prevent cross-contamination
- Patented safety heel inhibits blade tip from contacting handle
- Clean cool white LED illumination
- Low power use ensures long battery life
- Convenient for hospital and EMS use
- Compatible with all conventional ISO system handles
- Individually packaged; non-sterile

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5133-20	Premature	00	65 mm	12 mm	20
5-5133-00	Newborn	0	79 mm	12 mm	20
5-5133-01	Infant	1	105 mm	12 mm	20
5-5133-02	Child	2	157 mm	15 mm	20
5-5133-03	Medium Adult	3	200 mm	15 mm	20
5-5133-04	Large Adult	4	205 mm	17 mm	20

BLADE MEASUREMENTS



Conventional Laryngoscope

MACINTOSH AMERICAN

LATEX FREE | REUSABLE

SunMed


- Macintosh laryngoscope design is predominant choice among curved blades
- Smooth gentle spatula curve extends from the lock to the beak of the blade to indirectly raise the epiglottis effectively
- Blade design resembles that of a reverse Z-flange, facilitating greater access during intubation
- Frosted lamp and LED illumination available
- Frosted lamp provides defused illumination minus shadows
- LED lamp provides cool white illumination delivering 35,000 hours of use
- Constructed of 303/304 surgical stainless steel

FROSTED ILLUMINATION

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5052-00	Neonate	0	80 mm	15 mm	1
5-5052-01	Infant	1	87 mm	15 mm	1
5-5052-02	Child	2	108 mm	15 mm	1
5-5052-03	Medium Adult	3	128 mm	24 mm	1
5-5052-35	Ext. Med. Adult	3.5	144 mm	22 mm	1
5-5052-04	Large Adult	4	159 mm	22 mm	1
5-0234-01	Replacement Lamp for Sizes 0, 1, 2				10
5-0234-04	Replacement Lamp for Sizes 3, 3.5, 4				10

LED ILLUMINATION

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5852-00	Neonate	0	80 mm	15 mm	1
5-5852-01	Infant	1	87 mm	15 mm	1
5-5852-02	Child	2	108 mm	15 mm	1
5-5852-03	Medium Adult	3	128 mm	24 mm	1
5-5852-35	Ext. Med. Adult	3.5	144 mm	22 mm	1
5-5852-04	Large Adult	4	159 mm	22 mm	1
5-0234-21	Replacement Lamp for Sizes 0, 1, 2				1
5-0234-24	Replacement Lamp for Sizes 3, 3.5, 4				1



WATERPROOF

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5401-00	Neonate	0	80 mm	15 mm	1
5-5401-01	Infant	1	87 mm	15 mm	1
5-5401-02	Child	2	108 mm	15 mm	1
5-5401-03	Medium Adult	3	128 mm	24 mm	1
5-5401-35	Ext. Med. Adult	3.5	144 mm	22 mm	1
5-5401-04	Large Adult	4	159 mm	22 mm	1
5-5401-01	Replacement Lamp for Sizes 0, 1, 2				10
5-5401-04	Replacement Lamp for Sizes 3, 3.5, 4				10

Conventional Laryngoscope
SIKER MIRROR AMERICAN

LATEX FREE | REUSABLE




- Ideal for use on patients with anatomical variations which make intubation with conventional blades difficult or impossible
- The basic portions of the blade make an angle of 135° with a mirror located at flange angle facing spatula for better visualization
- Distal portion is 3"
- Mirror inverts the reflected image – clinician should exercise caution until familiar with its use
- Curved stylette is recommended for use with the Siker blade
- Immersing blade in warm water bath for 10 minutes prior to use will help prevent mirror fogging
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3075-04	Large Adult	4	158 mm	18 mm	1
5-0234-04	Replacement Lamp				10


Conventional Laryngoscope
IMPROVED VIEW AMERICAN

LATEX FREE | REUSABLE



- Channel extends from heel to mid-blade for improved visualization of the cords
- Channel design helps guide endotracheal tube insertion
- Clear reflector lamp delivers four times brighter illumination than frosted lamp
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5054-02	Child	2	100 mm	18 mm	1
5-5054-03	Medium Adult	3	130 mm	24 mm	1
5-5054-35	Ext. Med. Adult	3.5	144 mm	25 mm	1
5-5054-04	Ext. Lg. Adult	4	155 mm	26 mm	1
5-0234-14	Replacement Lamp				10



Conventional Laryngoscope

PRISMVIEW AMERICAN

LATEX FREE | REUSABLE

SunMed



- Designed with an optically polished prism with a 30° refraction without image inversion
- Ideal for use when the view of the vocal cords is impaired, especially the most anterior larynx, during nasal intubation, and for post operative examination of the larynx
- Help prevent condensation by immersing blade in warm water bath for 10 minutes prior to use
- Constructed of 303/304 surgical stainless steel


ITEM	DESCRIPTION	SIZE	A	B	PK
5-5029-02	Child	2	100 mm	17 mm	1
5-5029-03	Medium Adult	3	130 mm	24 mm	1
5-5029-35	Ext. Med. Adult	3.5	144 mm	25 mm	1
5-5029-04	Ext. Lg. Adult	4	155 mm	26 mm	1
5-0234-14	Replacement Lamp				5

Conventional Laryngoscope

FLANGELESS AMERICAN

LATEX FREE | REUSABLE

SunMed



- Bizzarri-Giuffrida
- Absent upper flange on curved blade is effective for use on patients with limited mouth opening, prominent incisors, receding mandible, short thick neck, or larynx in an extreme anterior anatomical position
- Design greatly reduces the potential trauma during laryngoscopy
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3072-03	Medium Adult	3	130 mm	8 mm	1
5-3072-04	Large Adult	4	158 mm	8 mm	1
5-0234-04	Replacement Lamp				10

Conventional Laryngoscope

BLECHMAN AMERICAN

LATEX FREE | REUSABLE

SunMed

- Angled tip elevates the epiglottis in supine patient
- Tongue is swept to the left allowing the blade to tip upward for effective viewing from the right side
- Clear lamp provides bright illumination
- Constructed of 303/304 surgical stainless steel



ITEM	DESCRIPTION	SIZE	A	B	PK
5-3084-03	Medium Adult	3	148 mm	14 mm	1
5-0234-04	Replacement Lamp				10

Conventional Laryngoscope

MACINTOSH ENGLISH

LATEX FREE | REUSABLE

SunMed

- Macintosh laryngoscope design is predominant choice among curved blades
- Flange extends all the way down to distal tip
- Soft matte finish virtually eliminates reflection and glare
- Clear lamp and LED illumination available
- Clear lamp provides bright illumination
- Cool white LED illumination delivers 35,000 hours of use; durable performance plus low current use increases battery life
- Constructed of 303/304 surgical stainless steel

CLEAR HALOGEN ILLUMINATION

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5051-00	Neonate	0	80 mm	17 mm	1
5-5051-01	Infant	1	92 mm	14 mm	1
5-5051-02	Child	2	100 mm	17 mm	1
5-5051-03	Medium Adult	3	130 mm	18 mm	1
5-5051-04	Large Adult	4	155 mm	18 mm	1
5-5051-05	X-Large Adult	5	175 mm	24 mm	1
5-0235-01	Replacement Lamp for Sizes 0, 1				10
5-0235-04	Replacement Lamp for Sizes 2, 3, 4, 5				10

LED ILLUMINATION

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5851-00	Neonate	0	80 mm	15 mm	1
5-5851-01	Infant	1	92 mm	14 mm	1
5-5851-02	Child	2	100 mm	17 mm	1
5-5851-03	Medium Adult	3	130 mm	18 mm	1
5-5851-04	Large Adult	4	155 mm	18 mm	1
5-5851-05	X-Large Adult	5	175 mm	24 mm	1
5-0234-21	Replacement Lamp for Sizes 0, 1				10
5-0234-24	Replacement Lamp for Sizes 2, 3, 4, 5				10




LEFT HANDED

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5053-03	Medium Adult Left	3	130 mm	18 mm	1
5-0235-04	Replacement Lamp				10



BLADE MEASUREMENTS



Conventional Laryngoscope

SUNFLEX™ ENGLISH

LATEX FREE | REUSABLE

SunMed



- English channel helps visualization the epiglottis
- Designed to enable precise control to elevate the epiglottis
- Non-pinching tip
- Satin finish virtually eliminates reflection and glare
- Reflector lamp provides extra bright illumination
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5082-02	Child	2	102 mm	20 mm	1
5-5082-03	Medium Adult	3	132 mm	25 mm	1
5-5082-04	Large Adult	4	155 mm	26 mm	1
5-0234-14	Replacement Lamp				5

Conventional Laryngoscope

PRISMVIEW ENGLISH

LATEX FREE | REUSABLE

SunMed




- Designed with an optically polished prism with a 30° refraction without image inversion
- Ideal for use when the view of the vocal cords is impaired, especially the most anterior larynx, during nasal intubation, and for post operative examination of the larynx
- Help prevent condensation by immersing blade in warm water bath for 10 minutes prior to use
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5028-02	Child	2	100 mm	17 mm	1
5-5028-03	Medium Adult	3	130 mm	24 mm	1
5-5028-35	Ext. Med. Adult	3.5	144 mm	25 mm	1
5-5028-04	Ext. Lg. Adult	4	155 mm	26 mm	1
5-0234-14	Replacement Lamp				5

Conventional Laryngoscope
REDUCED FLANGE E-MAC™

LATEX FREE | REUSABLE




- Reduced flange at rear of blade increases maneuverability and reduces pressure on teeth
- Channel extends from heel to mid-blade for improved visualization of the cords
- Channel design helps guide endotracheal tube insertion
- Clear reflector lamp delivers four times brighter illumination than frosted lamp
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3071-02	Child	2	100 mm	18 mm	1
5-3071-03	Medium Adult	3	130 mm	24 mm	1
5-3071-35	Ext. Med. Adult	3.5	144 mm	25 mm	1
5-3071-04	Ext. Lg. Adult	4	155 mm	26 mm	1
5-0234-14	Replacement Lamp				5


Conventional Laryngoscope
IMPROVED VIEW ENGLISH

LATEX FREE | REUSABLE



- Channel extends from heel to mid-blade for improved visualization of the cords
- Channel design helps guide endotracheal tube insertion
- Clear reflector lamp delivers four times brighter illumination than frosted lamp
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5046-02	Child	2	100 mm	18 mm	1
5-5046-03	Medium Adult	3	130 mm	24 mm	1
5-5046-35	Ext. Med. Adult	3.5	144 mm	25 mm	1
5-5046-04	Ext. Lg. Adult	4	155 mm	26 mm	1
5-0234-14	Replacement Lamp				5



Conventional Laryngoscope
POLIO ENGLISH

LATEX FREE | REUSABLE




- Modified Macintosh English profile blade for indirect laryngoscopy angled at 120°
- Blade is offset from handle at an obtuse angle to allow intubation of patients on respirators, in body jackets, as well as other difficult intubation scenarios
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3077-53	Medium Adult Angled at 120°	3	130 mm	24 mm	1
5-0235-04	Replacement Lamp				10

Conventional Laryngoscope
INTUBATION KIT

LATEX FREE



8-1050-85 Conventional Intubation Kit
Batteries, lubricating jelly, and syringe not shown.

ITEM	DESCRIPTION	PK
8-1050-85	Conventional Intubation Kit Includes items above + Conventional Laryngoscope Blades	1 Kit
8-1050-80	Nylon Case Only	1

Conventional Laryngoscope

MILLER AMERICAN

LATEX FREE | REUSABLE

SunMed

- Miller laryngoscope design is predominant choice among straight blades
- Reduced flange helps minimize trauma
- Distal end curve improves lifting of epiglottis
- Designed to facilitate greater exposure of larynx
- Blade sizes 2, 3, and 4 have lamps on right side
- Frosted lamp and LED illumination available
- Frosted lamp provides defused illumination without shadows
- Cool white LED illumination delivers 35,000 hours of use; durable performance plus low current use increases battery life
- Constructed of 303/304 surgical stainless steel




FROSTED HALOGEN ILLUMINATION

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5062-20	Premature	00	65 mm	10 mm	1
5-5062-00	Neonate	0	75 mm	10 mm	1
5-5062-01	Infant	1	102 mm	10 mm	1
5-5062-15	Small Child	1.5	130 mm	12 mm	1
5-5062-02	Child	2	155 mm	13 mm	1
5-5062-03	Medium Adult	3	195 mm	13 mm	1
5-5062-04	Large Adult	4	205 mm	14 mm	1
5-0234-01	Replacement Lamp for Sizes 00, 0, 1, 1.5, 2				10
5-0234-04	Replacement Lamp for Sizes 3, 4				10

LED ILLUMINATION

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5862-20	Premature	00	65 mm	10 mm	1
5-5862-00	Neonate	0	75 mm	10 mm	1
5-5862-01	Infant	1	102 mm	10 mm	1
5-5862-15	Small Child	1.5	130 mm	12 mm	1
5-5862-02	Child	2	155 mm	13 mm	1
5-5862-03	Medium Adult	3	195 mm	13 mm	1
5-5862-04	Large Adult	4	205 mm	14 mm	1
5-0234-21	Replacement Lamp for Sizes 00, 0, 1, 1.5, 2				1
5-0234-24	Replacement Lamp for Sizes 3, 4				1



WATERPROOF

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5403-20	Premature	00	65 mm	10 mm	1
5-5403-00	Neonate	0	75 mm	10 mm	1
5-5403-01	Infant	1	102 mm	10 mm	1
5-5403-15	Small Child	1.5	130 mm	12 mm	1
5-5403-02	Child	2	155 mm	13 mm	1
5-5403-03	Medium Adult	3	195 mm	13 mm	1
5-5403-04	Large Adult	4	205 mm	14 mm	1
5-0234-01	Replacement Lamp for Sizes 00, 0, 1, 1.5, 2				10
5-0234-04	Replacement Lamp for Sizes 3, 4				10

Conventional Laryngoscope
SNOW AMERICAN

LATEX FREE | REUSABLE



- Combined Miller spatula and the Wis-Foregger flange design can reduce trauma during intubation, especially for use with a cuffed tube
- Curves back from the distal end for a convenient visual field
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3069-03	Medium Adult	3	162 mm	16 mm	1
5-0234-04	Replacement Lamp				10


Conventional Laryngoscope
CRANWALL AMERICAN

LATEX FREE | REUSABLE




- Reduced flange allows easier insertion through a restricted opening, such as arthritic jaw or facial burns
- Design helps decrease potential damage to upper teeth
- Flange removed directly below the light carrier
- Can easily be used as a left-handed straight blade
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5063-00	Neonate	0	75 mm	6 mm	1
5-5063-01	Infant	1	100 mm	6 mm	1
5-5063-02	Child	2	157 mm	8 mm	1
5-5063-03	Medium Adult	3	195 mm	9 mm	1
5-0234-01	Replacement Lamp for Sizes	0, 1, 2			10
5-0234-04	Replacement Lamp for Sizes	3			10



Conventional Laryngoscope
PHILLIPS AMERICAN

LATEX FREE | REUSABLE




- Combines the straight Jackson blade design with the curved distal tip Miller design to allow greater visibility and a 'direct line' approach to the trachea
- Low profile flange reduces the risks of oral damage
- Unique lamp mounting enables deep illumination both downward and inward
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3073-01	Child	1	113 mm	13 mm	1
5-3073-02	Adult	2	166 mm	18 mm	1
5-0234-01	Replacement Lamp				10


Conventional Laryngoscope
WISCONSIN AMERICAN

LATEX FREE | REUSABLE



- Straight spatula design with flange that expands slightly toward the distal portion of the blade increases visual field helping to prevent potential trauma during intubation
- Distal portion of the blade is wider and formed slightly to the right, making it well-adapted to lifting the epiglottis
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5057-02	Child	2	135 mm	17 mm	1
5-5057-03	Medium Adult	3	162 mm	18 mm	1
5-5057-04	Large Adult	4	199 mm	21 mm	1
5-0234-04	Replacement Lamp				10



Conventional Laryngoscope

WIS-HIPPLE AMERICAN

LATEX FREE | REUSABLE

SunMed



- Wisconsin-modified blade with a straight spatula
- Flange tends to be straighter than a Wisconsin blade and runs parallel with the spatula
- Designed primarily for use in infants
- Wider, flat tip is well-adapted to lifting the epiglottis
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5057-20	Premature	00	65 mm	11 mm	1
5-5057-00	Newborn	0	76 mm	12 mm	1
5-5057-01	Infant	1	102 mm	10 mm	1
5-5057-15	Small Child	1.5	115 mm	14 mm	1
5-0234-01	Replacement Lamp				10

Conventional Laryngoscope

MILLER PORT AMERICAN

LATEX FREE | REUSABLE

SunMed



- Specifically designed to perform laryngoscopy and intubation with reduced risk of hypoxia
- Integrated tube permits delivery of oxygen or other gas mixtures during intubation
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3078-00	Newborn	0	75 mm	10 mm	1
5-3078-01	Infant	1	102 mm	10 mm	1
5-0234-04	Replacement Lamp				10

Conventional Laryngoscope
SCHAPIRA AMERICAN

LATEX FREE | REUSABLE



- Absent flange helps reduces the risk of trauma to the patient while providing excellent visualization
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3080-04	Large Adult	4	177 mm	7 mm	1
5-0234-04	Replacement Lamp				10


Conventional Laryngoscope
WHITEHEAD AMERICAN

LATEX FREE | REUSABLE



- Reduced flange design increases visualization and facilitates intubation
- Helps prevent potential damage to the incisors
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel


ITEM	DESCRIPTION	SIZE	A	B	PK
5-3065-03	Medium Adult	3	162 mm	13 mm	1
5-0234-04	Replacement Lamp				10



Conventional Laryngoscope MILLER ENGLISH

LATEX FREE | REUSABLE

SunMed



- Miller laryngoscope design is predominant choice among straight blades
- Satin finish virtually eliminates reflection and glare
- Halogen and LED illumination available
- Clear halogen provides bright illumination
- Cool white LED illumination delivers 35,000 hours of use; durable performance plus low current use increases battery life
- Constructed of 303/304 surgical stainless steel

Conventional Laryngoscope OXFORD ENGLISH

LATEX FREE | REUSABLE

SunMed



- Suitable for premature infants, babies and young children up to the age of four
- Overhang on the open side prevents lips from obscuring vision
- Broad, flat surface aids pediatric patients with extreme cleft palates
- Clear lamp provides bright illumination
- Constructed of 303/304 surgical stainless steel

CLEAR HALOGEN ILLUMINATION

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5061-00	Neonate	0	77 mm	9 mm	1
5-5061-01	Infant	1	102 mm	9 mm	1
5-5061-02	Child	2	153 mm	11 mm	1
5-5061-03	Medium Adult	3	195 mm	11 mm	1
5-5061-04	Large Adult	4	205 mm	16 mm	1
5-0235-01	Replacement Lamp	0, 1			10
5-0235-04	Replacement Lamp	2, 3, 4			10

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5078-01	Infant	1	102 mm	12 mm	1
5-0235-01	Replacement Lamp				10

LED ILLUMINATION

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5861-00	Neonate	0	77 mm	9 mm	1
5-5861-01	Infant	1	102 mm	9 mm	1
5-5861-02	Child	2	153 mm	11 mm	1
5-5861-03	Medium Adult	3	195 mm	11 mm	1
5-5861-04	Large Adult	4	205 mm	16 mm	1
5-0234-21	Replacement Lamp	0, 1			1
5-0234-14	Replacement Lamp	2, 3, 4			1

SunMed

FOR COMPLETE
Customer Service
Technical Information
Care/Cleaning Instructions
Product Literature

CALL
800-433-2797

VISIT
www.Sun-Med.com

EMAIL
info@Sun-Med.com

Conventional Laryngoscope
GUEDEL ENGLISH

LATEX FREE | REUSABLE



- 28° angled blade helps promote lifting without using the teeth as a fulcrum
- Designed specifically for use with cuffed endotracheal tubes
- Distal tip is angled slightly more than the Flagg to assist in epiglottis compression
- Frosted lamp provides defused illumination without shadows
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-3055-01	Infant	1	97 mm	13 mm	1
5-3055-02	Child	2	122 mm	13 mm	1
5-3055-03	Medium Adult	3	137 mm	16 mm	1
5-3055-04	Large Adult	4	162 mm	17 mm	1
5-0234-01	Replacement Lamp	1, 2			10
5-0234-04	Replacement Lamp	3, 4			10


Conventional Laryngoscope
ROBERTSHAW ENGLISH

LATEX FREE | REUSABLE



- Designed specifically for infants and children
- Spatula is gently curved over the distal third and can indirectly lift the epiglottis similar to a Macintosh blade
- Permits binocular vision allowing better judgment of depth and consequently less risk of trauma
- Clear lamp provides bright illumination
- Constructed of 303/304 surgical stainless steel


ITEM	DESCRIPTION	SIZE	A	B	PK
5-5077-00	Newborn	0	94 mm	10 mm	1
5-5077-01	Infant	1	112 mm	12 mm	1
5-0235-01	Replacement Lamp				1



Conventional Laryngoscope
SEWARD ENGLISH

LATEX FREE | REUSABLE

SunMed



- Originally produced as an infant blade for use in asphyxia neonatorum, it is now available to fit older children
- Clear lamp provides bright illumination
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5079-01	Infant	1	105 mm	7 mm	1
5-5079-02	Child	2	130 mm	13 mm	1
5-0235-01	Replacement Lamp				10

Conventional Laryngoscope
SOPER ENGLISH

LATEX FREE | REUSABLE

SunMed



- Broad flat blade makes it easier to restrict premature and newborn infant tongue movement
- Clear lamp provides bright illumination
- Constructed of 303/304 surgical stainless steel

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5080-01	Infant	1	100 mm	10 mm	1
5-0235-01	Replacement Lamp				10

Conventional Laryngoscope

LARYNGOSCOPE KITS

LATEX FREE | REUSABLE



- Available in Macintosh and Miller, American and English profiles
- Blades constructed of 303/304 surgical stainless steel
- High impact plastic case - for ease of transport
- Include most popular laryngoscope blades and replacement lamps
- Include chrome plated handle



MACINTOSH AMERICAN PROFILE

ITEM	KIT INCLUDES	PK
5-5052-44	Medium Handle (5-0237-03) Macintosh Blades: 1 (5-5052-01), 2 (5-5052-02), 3 (5-5052-03), 4 (5-5052-04) 2 Replacement Lamps (5-0234-01, 5-0234-04) Case (8-1050-73)	1 Kit



MACINTOSH ENGLISH PROFILE

ITEM	KIT INCLUDES	PK
5-5051-44	Medium Handle (5-0237-03) Macintosh Blades: 1 (5-5051-01), 2 (5-5051-02), 3 (5-5051-03), 4 (5-5051-04) 2 Replacement Lamps (5-0235-01, 5-0235-04) Case (8-1050-73)	1 Kit



MILLER ENGLISH PROFILE

ITEM	KIT INCLUDES	PK
5-5061-45	Medium Handle (5-0237-03) Miller Blades: 0 (5-5061-00), 1 (5-5061-01), 2 (5-5061-02), 3 (5-5061-03), 4 (5-5061-04) 2 Replacement Lamps (5-0235-01, 5-0235-04) Case (8-1050-71)	1 Kit



MILLER AMERICAN PROFILE

ITEM	KIT INCLUDES	PK
5-5062-45	Medium Handle (5-0237-03) Miller Blades: 0 (5-5062-00), 1 (5-5062-01), 2 (5-5062-02), 3 (5-5062-03), 4 (5-5062-04) 2 Replacement Lamps (5-0234-01, 5-0234-04) Case (8-1050-71)	1 Kit



MACINTOSH/MILLER AMERICAN PROFILE

ITEM	KIT INCLUDES	PK
5-5062-47	Medium Handle (5-0237-03) Macintosh Blades: 2 (5-5052-02), 3 (5-5052-03), 4 (5-50522-04) Miller Blades: 0 (5-5062-00), 1 (5-5062-01), 2 (5-5062-02), 3 (5-5062-03) 2 Replacement Lamps (5-0234-01, 5-0234-04) Case (8-1050-77)	1 Kit



MACINTOSH/MILLER ENGLISH PROFILE

ITEM	KIT INCLUDES	PK
5-5061-47	Medium Handle (5-0237-03) Macintosh Blades: 2 (5-5051-02), 3 (5-5051-03), 4 (5-50521-04) Miller Blades: 0 (5-5061-00), 1 (5-5061-01), 2 (5-5061-02), 3 (5-5061-03) 2 Replacement Lamps (5-0235-01, 5-0235-04) Case (8-1050-77)	1 Kit



SUN-FLEX

ITEM	KIT INCLUDES	PK
5-5082-42	Medium Handle (5-0237-03) SunFlex Blades: 3 (5-5082-03), 4 (5-5082-04) Replacement Lamp (5-0234-14) Case (8-1050-72)	1 Kit
8-1050-72	SunFlex Blade Case Only	1

Conventional Laryngoscope
CHROME PLATED

LATEX FREE | REUSABLE



- Chrome plated brass
- Batteries are not included

ITEM	DESCRIPTION	SIZE	PK	BATTERIES
5-0237-01	Chrome Plated Brass	Penlite	1	2 "AA"
5-0237-02	Chrome Plated Brass	Stubby	1	2 "AA"
5-0237-03	Chrome Plated Brass	Medium	1	2 "C"
5-0237-05	Chrome Plated Brass	Large	1	2 "D"
5-0237-08	Chrome Plated Brass	Mini	1	Lithium 123

Conventional Laryngoscope
STAINLESS STEEL WATERPROOF

LATEX FREE | REUSABLE



- Stainless steel waterproof handles may be immersed in disinfectant without removing batteries
- Batteries are not included

ITEM	DESCRIPTION	SIZE	PK	BATTERIES
5-0237-91	Waterproof Handle	Penlite	1	2 "AA"
5-0237-92	Waterproof Handle	Stubby	1	2 "AA"
5-0237-93	Waterproof Handle	Medium	1	2" C"
5-0237-96	Waterproof Handle	Micro	1	Lithium 123

Conventional Laryngoscope
HOWLAND LOCK

LATEX FREE | REUSABLE



- Natural lifting simplifies laryngeal exposure, making intubation possible even in the most difficult cases
- Built-in leverage decreases prying and reduces possibility of broken teeth
- Compatible with all hook-on conventional laryngoscope handles and blades
- Can often make an otherwise certain tracheotomy unnecessary

ITEM	DESCRIPTION	PK
5-0236-00	Conventional Howland Lock	1

Conventional Laryngoscope
DISPOSABLE LARYNGOSCOPE KIT

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE



- Ideal for ambulances, paramedics and first responders
- Variety of Macintosh and Miller blades
- Available with medium conventional handles
- Components arrive in original packaging inside sturdy canvas pouch

ITEM	KIT INCLUDES	PK
5-5133-57	Medium Handle (5-0237-03) Macintosh Blades: 2 (5-5132-02), 3 (5-5132-03), 4 (5-5132-04) Miller Blades: 0 (5-5133-00), 1(5-5133-01), 2 (5-5133-02), 3 (5-5133-03) Canvas Pouch	1

Conventional Laryngoscope
REPLACEMENT BLADE LAMPS

LATEX FREE | REUSABLE



- For use in conventional blades
- Compatible with all conventional blades
- Meets stringent manufacturing controls
- Frosted lamps provide defused illumination without shadows
- Clear lamps provide bright illumination
- Contact Customer Care for correct lamp selection

ITEM	SIZE	DESCRIPTION	PK
5-0234-01	Small	Frosted Lamp	10
5-0234-04	Large	Frosted Lamp	10
5-0234-14	Large	Reflector Lamp	5
5-0234-21	Small	LED Lamp	1
5-0234-24	Large	LED Lamp	1
5-0235-01	Small	Clear Lamp	10
5-0235-04	Large	Clear Lamp	10



Conventional Laryngoscope
INTUBATION KIT

LATEX FREE

SunMed



Complete intubation kits available with conventional,
fiber optic or fiber optic disposable laryngoscope blades

Each kit includes:

- 8 Oropharyngeal color coded Guedel airways (sizes 40 mm-100 mm)
- 1 Penlite chrome plated laryngoscope handle
- 1 Medium chrome plated laryngoscope handle
- 5 Miller laryngoscope blades (sizes 0-4)
- 4 Macintosh laryngoscope blades (sizes 1-4)
- 4 Extra laryngoscope lamps (2 small, 2 large)
- 1 Adult Magill forceps
- 1 Pediatric Magill forceps
- 2 “C” batteries
- 2 “AA” batteries
- 4 Packets of sterile lubricating jelly
- 2 Endotracheal uncuffed tubes (sizes 2.5, 4 mm)
- 8 Endotracheal cuffed tubes (sizes 5-10 mm)
- 4 Nasopharyngeal airways (sizes 20, 24, 28, 32 Fr)
- 1 10 mL syringe
- 1 Endotracheal tube introducer (bougie) with coude tip
- 1 Adult stylette
- 1 Nylon airway management case


ITEM	DESCRIPTION	PK
8-1050-95	Fiber Optic Intubation Kit Includes items above + Reusable FO Laryngoscope Blades	1
8-1050-85	Conventional Intubation Kit Includes items above + Conventional Laryngoscope Blades	1
8-1050-94	Single-Use Fiber Optic Intubation Kit Includes items above + Disposable FO Single Use Laryngoscope Blades	1
8-1050-80	Nylon Case Only	1

Inflatable Anesthesia Mask

VENT MASK™

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ventlab



- Superior ergonomic cushion design for an optimal seal with minimum pressure
- Flexible body allows gentle, deep seating of cushion to face
- Inflation valve adjusts cuff pressure for optimal seal
- Ample nose bridge room for proper fit
- Centered connector port for secure gripping and easy connect/disconnect
- Removable hook ring enables manual or head strap application
- Premium soft cushioned mask
- Pliable, ultra-thin wall
- Low-pressure soft seal with superior seating characteristics

ITEM	DESCRIPTION	SIZE	PK
VR1100	Vent Mask	Small Infant	20, 40
VR2100	Vent Mask	Infant	20, 40
VR3100	Vent Mask	Small Child	20, 40
VR4100	Vent Mask	Child	20, 40
VR5000	Vent Mask	Wide Adult	20, 50, 100
VR5100	Vent Mask	Adult	20, 50, 100
VR6100	Vent Mask	Large Adult	20, 50, 100

Inflatable Anesthesia Mask

SOFTAIRE™

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ventlab



- Non-conductive anesthesia mask
- Inflatable, soft cushioned mask
- Superior patient comfort
- Effective seal and seating

ITEM	DESCRIPTION	SIZE	PK
VR1700	Soft Aire Mask	Small Infant	20, 40
VR2700	Soft Aire Mask	Infant	20, 40
VR3700	Soft Aire Mask	Small Child	20, 40
VR4700	Soft Aire Mask	Child	20, 40
VR5700	Soft Aire Mask	Adult	20, 50, 100
VR6700	Soft Aire Mask	Wide Adult	20, 50, 100

Inflatable Mask

FLEXIBLE

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE



- Inflatable cuff helps ensure patient comfort
- Reduces wasted anesthetic gas pollution
- Transparent housing for immediate detection of regurgitation or silent aspiration
- Removable color-coded hook ring

ITEM	DESCRIPTION	COLOR	SIZE	PK
4-0110-00	Flexible Face Mask	Gray	Neonate	30
4-0110-01	Flexible Face Mask	Light Blue	Infant	30
4-0110-02	Flexible Face Mask	Yellow	Child	30
4-0110-03	Flexible Face Mask	Green	Small Adult	30
4-0110-04	Flexible Face Mask	Red	Medium Adult	30
4-0110-05	Flexible Face Mask	Dark Blue	Large Adult	30

Accessory

MASK HARNESS

LATEX FREE | REUSABLE



- Molded, not stamped
- Reinforced holes on adult size for extra durability
- Neoprene
- Center vents for comfort and fit at back of head
- Individually wrapped

ITEM	DESCRIPTION	SIZE	PK
4-1054-21	Infant Mask Harness	13" x 2.5"	1
4-1054-22	Child Mask Harness	22" x 4"	1
4-1054-23	Adult Mask Harness	27" x 5.5"	1

Non-Inflatable Anesthesia Mask

VENT MASK II™

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE



- Non-conductive anesthesia mask
- Air-filled, soft cushioned mask
- Pliable, ultra-thin wall
- Low-pressure soft seal with superior seating characteristics

ITEM	DESCRIPTION	SIZE	PK
VR0040	Vent Mask II	Neonate	20, 40
VR0400	Vent Mask II	Small Neonate	20, 40
VR1400	Vent Mask II	Small Infant	20, 40
VR2400	Vent Mask II	Infant	20, 40
VR3400	Vent Mask II	Small Child	20, 40
VR4400	Vent Mask II	Child	20, 40
VR5400	Vent Mask II	Adult	20, 50, 100
VR6400	Vent Mask II	Large Adult	20, 50, 100

Non-Inflatable Anesthesia Mask

SOFT SEAL™

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE



- Non-conductive anesthesia mask
- Non-inflatable mask
- Pliable cushion
- Economical anesthesia mask

ITEM	DESCRIPTION	SIZE	PK
VR1600	Soft Seal Mask	Neonate	20
VR2600	Soft Seal Mask	Infant	20
VR3600	Soft Seal Mask	Child	20

Rapid Infusion

INFU-SURG® PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox




INFU-SURG®

- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Pressure relief valve prevents over inflation (330 mm Hg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Color-coded gauge makes for accurate pressure monitoring (150-300 mm Hg)
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

ITEM	DESCRIPTION	VOLUME	PK
4005H	Infu-Surg® with Hook	500 mL	5
4010H	Infu-Surg® with Hook	1000 mL	5
4030	Infu-Surg®	3000 mL	5
4030H	Infu-Surg® with Hanger Loop	3000 mL	5

Rapid Infusion

INFU-SURG® CLARITY™ PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox



Clear panel enables full visibility of solution bag to monitor fluid level and easily scan patient barcode



INFU-SURG® CLARITY™

- Transparent panel for clear visibility of solution bag
- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Dedicated hook enables quick change of solution bag while providing proper fluid bag positioning in pocket
- Pressure relief valve prevents over-inflation (330 mm Hg)
- Color-coded gauge provides accurate pressure monitoring (150-300 mm Hg)
- Single-handed inflation and deflation with oval-shaped bulb makes it easy to use and requires minimal training
- Welded panel and hook construction are more durable than sewn designs to withstand repeated use
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

ITEM	DESCRIPTION	VOLUME	PK
4050	Infu-Surg® Clarity	500 mL	5
4100	Infu-Surg® Clarity	1000 mL	25

Rapid Infusion

VENTLAB PRESSURE INFUSION

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ventlab



Back mesh panel allows easy viewing of fluid bag.



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook integrated with bag loop to enable easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use.



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

- Pressure relief valve helps prevent over-inflation (330 mm Hg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Color-coded dial gauge makes for accurate pressure monitoring (0-300 mm Hg)
- Three-way stopcock ensures precise control of pressure
- Economical and dependable pressure infusion choice

ITEM	DESCRIPTION	SIZE	PK
VC405	Pressure Infusion	500 mL	25
VC410	Pressure Infusion	1000 mL	25
VC430	Pressure Infusion	3000 mL	25

Patient Safety

IGUARD EYE PROTECTOR

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE




- Rigid, clear protective patient eye cover for use in OR, ER and during general procedures
- Non-allergenic self-adhesive foam cushion for fast and accurate application
- Thick foam for patient comfort even in prone position
- Improved one piece backing with “courtesy tabs”

STERILE


ITEM	DESCRIPTION	PK
9-0210-00	Adult Eye Protector	25
9-0210-21	Pediatric Eye Protector	25
9-0210-20	Double Foam Eye Protector	25

NON-STERILE


ITEM	DESCRIPTION	PK
9-0210-00I	Adult Eye Protector	25
9-0210-21I	Pediatric Eye Protector	25
9-0210-20I	Double Foam Eye Protector	25



9-0210-00
Standard Foam



9-0210-20
Double Foam



Non-allergenic self-adhesive foam for easy application and removal



FOR COMPLETE
Customer Service
Technical Information
Care/Cleaning Instructions
Product Literature

CALL
800-433-2797

VISIT
www.Sun-Med.com

EMAIL
info@Sun-Med.com

Peripheral Nerve Stimulator

SUNSTIM PERIPHERAL NERVE STIMULATOR

LATEX FREE



SunStim™

8-1053-60 Features:

- Constant current stimulation for monitoring neuromuscular block
- Output current is fully adjustable from 0-70 mA

8-1053-60 Switch Functions:

- Stand-by position: no stimulus pulses are generated
- Twitch: 1 pulse per second
- Tetanus: 100 Hz (100 pulses/s)
- Built-in train of four: 4 pulses in a period of 2 seconds



SunStim™ Plus

8-1053-62 Features:

- Delivers sufficient output to ensure supramaximal stimulation; output current is fully adjustable from 0-70 mA

8-1053-62 Switch Functions:

- Double burst: 2 60-ms bursts of 50 Hz separated by .75 s
- Twitch: 1 pulse per second
- Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
- Built-in train of four: 4 pulses in a period of 2 seconds
- Stimulus pulse: square wave monophasic pulses (200-μs duration)



All units include instruction manual, ball electrodes, battery (9V alkaline) and 6' lead wires with alligator clips

PERIPHERAL NERVE STIMULATOR

ITEM	DESCRIPTION	PK
8-1053-60	SunStim™ Peripheral Nerve Stimulator	1
8-1053-62	SunStim™ Plus Peripheral Nerve Stimulator	1
8-1053-63	SunStim™ Pro Peripheral Nerve Stimulator	1

- Delivers reliable, adjustable current in one-second twitch pulses
- Features built-in train of four with four pulses per two seconds
- Stimulator output current is fully adjustable from 0-70 mA
- Visual illumination for pulse activity and low battery life
- Lightweight and pocket-size for easy use and carry

SunStim™ Pro

8-1053-63 Features:

- Easy-to-read digital screen display for DBS, Tetanus, TOF and Twitch functionality
- Ergonomic and compact shape; fits better in hand and more comfortable to use
- Includes removable gray silicone protective wrap

8-1053-63 Switch Functions:

- Double burst: 2 60-ms bursts of 50 Hz separated by .75s
- Twitch: 1 pulse per second
- Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
- Built-in train of four: 4 pulses in a period of 2 seconds
- Stimulus pulse: square wave monophasic pulses (200-μs duration)



Removable Silicone Pro Wrap helps protect device, while providing a secure grip




8-1053-76



8-1053-78



8-1053-79



8-1053-85
8-1053-95



8-1053-87

SILICONE PRO WRAP (FOR USE ONLY WITH 8-1053-63)

ITEM	DESCRIPTION	COLOR	PK
8-1053-GR	Silicone Pro Wrap	Grey	1
8-1053-BL	Silicone Pro Wrap	Blue	1
8-1053-PL	Silicone Pro Wrap	Purple	1
8-1053-FS	Silicone Pro Wrap	Fuchsia	1
8-1053-PK	Silicone Pro Wrap	Pink	1
8-1053-RD	Silicone Pro Wrap	Red	1
8-1053-OR	Silicone Pro Wrap	Orange	1
8-1053-GN	Silicone Pro Wrap	Green	1
8-1053-TL	Silicone Pro Wrap	Teal	1

ACCESSORIES

ITEM	DESCRIPTION	PK
8-1053-76	Ball Electrode and Hollow Pin	1
8-1053-78	6' Lead, Snap and Shielded Pin	1
8-1053-79	6' Lead, Clip and Shielded Pin	1
8-1053-85	Pole Mount Bracket for 8-1053-60/-62	1
8-1053-95	Pole Mount Bracket for 8-1053-63	1
8-1053-87	Zippered Pouch	1

Monoscope

COMFORT FIT SOFTEAR™

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

SunMed

- Designed specifically to address sound transmission challenges for the anesthesiologist due to the irregular and variable shape of the ear canal
- Once compressed foam earpiece is inserted, it expands slowly for a perfect, low pressure seal with the canal wall
- Sealing action enables the accurate transmission of low volume sounds
- Low pressure fit against ear canal provides maximum comfort, thus promoting more vigilant use
- Connects to esophageal stethoscopes and Wenger (precordial) chest pieces
- 4' Tubing

ITEM	DESCRIPTION	PK
6-1030-25	Comfort FIT SoftEar™ Monoscope	25

Precordial Stethoscope

CHEST PIECES

LATEX FREE

SunMed

- Durable stainless steel construction – no chrome plating to flake or peel
- Pediatric chest pieces are hollow and lightweight, with rounded, smooth chest surface, making them ideal for use with newborns and infants
- Solid Wenger chest pieces feature a relatively deep sound chamber
- All chest pieces have a luer slip male connector
- Reusable

ITEM	DESCRIPTION	SIZE	PK
6-0025-00	Newborn Chest Piece	0.875" Diameter	1
6-0025-01	Infant Chest Piece	1.125" Diameter	1
6-0025-02	Child Chest Piece	1.125" Diameter	1
6-0025-03	Adult Wenger Chest Pc	1.44" Diameter	1

Precordial Stethoscope

PRECORDIAL SET

LATEX FREE | SINGLE PATIENT USE

SunMed

- 1 each adult, child, infant, and newborn chest piece
- 8 Double adhesive disks (also available separately)
- 1 Comfort FIT SoftEar™ disposable monoscope



ITEM	DESCRIPTION	PK
6-0025-05	Precordial Stethoscope Set	1
6-0025-25	Double Adhesive Disks	100

Stethoscope SPRAGUE

LATEX FREE

SunMed



- Three different size bells eliminates the need for multiple stethoscopes for infants, children and adults
- Various heads for acoustical responses in low, medium and high pitch sounds
- Heavy wall 22" tubing and choice of ear tips help reduce sound leakage and extraneous noise
- Complete stethoscope kit includes 1 pair black plastic ear tips, 3 open bells (infant, child, adult), 1 pair soft rubber (mushroom type) ear tips, adult and pediatric diaphragms
- 29" Overall length

Stethoscope DUAL HEAD

LATEX FREE

SunMed



- Molded one piece Y tubing
- High sensitivity - lightweight diaphragm
- Non-chill ring
- Anodized aluminum chest-piece
- Chrome plated brass binaurals
- White plastic ear tips
- Superior acoustics
- 31" Overall length

ITEM	DESCRIPTION	PK
6-0028-00	Sprague Rappaport Stethoscope	1

ITEM	DESCRIPTION	COLOR	PK
6-0027-40	Dual Head Stethoscope	Black	1
6-0027-41	Dual Head Stethoscope	Red	1
6-0027-42	Dual Head Stethoscope	Grey	1
6-0027-43	Dual Head Stethoscope	Green	1

Oximeter

FINGERTIP PULSE OXIMETER

LATEX FREE

SunMed

- Accurately displays patient O₂ saturation and pulse rate using high readability OLED display technology
- SPO₂ measurement from 70-99%; accuracy 80-90% +/- 2%
- Pulse rate measurement 30-225 BPM, +/- 2BPM, 100-235 BPM +/- 2%
- AAA battery life of up to 30 hours; auto shut-off feature, low power indicator, 10 brightness scales
- Convenient pocket size: 2" x 1.25" x 1.375"; wt. 1.2 oz.

ITEM	DESCRIPTION	PK
9-6001-00	Fingertip Pulse Oximeter	1



Includes lanyard, protective cover and carrying pouch