



EMERGENCY




Oropharyngeal Airway

BERMAN

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Firm, polyethylene plastic
- Vented septum
- Built-in bite block
- Sets are ideal for crash carts
- Individually wrapped

INDIVIDUAL

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1508-40	Newborn	Pink	4 cm/ 40 mm	10
1-1508-50	Infant	Blue	5 cm/ 50 mm	10
1-1508-60	Small Child	Black	6 cm/ 60 mm	10
1-1508-70	Child	White	7 cm/ 70 mm	10
1-1508-80	Small Adult	Green	8 cm/ 80 mm	10
1-1508-90	Medium Adult	Yellow	9 cm/ 90 mm	10
1-1508-99	Large Adult	Red	10 cm/ 100 mm	10
1-1508-11	X-Large Adult	Orange	11 cm/ 110 mm	10

1-1506-40	Newborn	Pink	4 cm/ 40 mm	50
1-1506-50	Infant	Blue	5 cm/ 50 mm	50
1-1506-60	Small Child	Black	6 cm/ 60 mm	50
1-1506-70	Child	White	7 cm/ 70 mm	50
1-1506-80	Small Adult	Green	8 cm/ 80 mm	50
1-1506-90	Medium Adult	Yellow	9 cm/ 90 mm	50
1-1506-99	Large Adult	Red	10 cm/ 100 mm	50
1-1506-11	X-Large Adult	Orange	11 cm/ 110 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-5008-00	40, 50, 60, 70, 80, 90, 100, 110 mm (1 of each size per set; 10 sets per pack)	10 Set

Oropharyngeal Airway

TRADITIONAL GUEDEL

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Firm, polyethylene plastic
- One piece construction
- Built-in bite block
- Sets are ideal for crash carts
- Bulk packaged

INDIVIDUAL

ITEM	DESCRIPTION	COLOR	SIZE	PK
1-1501-60	Small Child	Black	6 cm/ 60 mm	10
1-1501-70	Child	White	7 cm/ 70 mm	10
1-1501-80	Small Adult	Green	8 cm/ 80 mm	10
1-1501-90	Medium Adult	Yellow	9 cm/ 90 mm	10

1-1500-60	Small Child	Black	6 cm/ 60 mm	50
1-1500-70	Child	White	7 cm/ 70 mm	50
1-1500-80	Small Adult	Green	8 cm/ 80 mm	50
1-1500-90	Medium Adult	Yellow	9 cm/ 90 mm	50

SETS

ITEM	SIZES INCLUDED	PK
1-1501-00	60, 70, 80, 90 mm (1 of each size per set/pack)	1 Set

Nasopharyngeal Airway

LATEX FREE ROBERTAZZI

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Blunt tip
- Resists collapse or kinking
- Packaged with large size and specification notation for accurate selection
- Sets are ideal for crash carts
- Individually sterile packaged

INDIVIDUAL

ITEM	SIZE	ID	LENGTH	PK
1-5075-20	20 Fr	3.9 mm	114 mm	10
1-5075-22	22 Fr	4.8 mm	114 mm	10
1-5075-24	24 Fr	5.3 mm	120 mm	10
1-5075-26	26 Fr	6.3 mm	121 mm	10
1-5075-28	28 Fr	6.8 mm	132 mm	10
1-5075-30	30 Fr	7.2 mm	139 mm	10
1-5075-32	32 Fr	7.6 mm	146 mm	10
1-5075-34	34 Fr	8.5 mm	157 mm	10
1-5075-36	36 Fr	9.0 mm	170 mm	10

SETS

ITEM	SIZES INCLUDED (1 of each size per set/pack)	PK
1-5075-04	20, 24, 28, 32 Fr + sterile lubricating jelly	1 Set
1-5075-05	28, 30, 32, 34, 36 Fr + sterile lubricating jelly	1 Set
1-5075-06	22, 24, 26, 28, 30, 32 Fr + sterile lubricating jelly	1 Set
1-5075-09	20, 22, 24, 26, 28, 30, 32, 34, 36 Fr + sterile lubricating jelly	1 Set

Nasopharyngeal Airway

ADJUSTABLE FLANGE AIRWAY

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Soft, adjustable flange
- Blunt tip
- Resists collapse or kinking
- Individually sterile packaged

INDIVIDUAL

ITEM	SIZE	ID	LENGTH	PK
1-5072-12	12 Fr	3.0 mm	95 mm	10
1-5072-14	14 Fr	3.5 mm	95 mm	10
1-5072-16	16 Fr	4.0 mm	95 mm	10
1-5072-18	18 Fr	4.5 mm	95 mm	10
1-5072-20	20 Fr	5.0 mm	115 mm	10
1-5072-22	22 Fr	5.5 mm	120 mm	10
1-5072-24	24 Fr	6.0 mm	120 mm	10
1-5072-26	26 Fr	6.5 mm	130 mm	10
1-5072-28	28 Fr	7.0 mm	170 mm	10
1-5072-30	30 Fr	7.5 mm	170 mm	10
1-5072-32	32 Fr	8.0 mm	170 mm	10
1-5072-34	34 Fr	8.5 mm	170 mm	10
1-5072-36	36 Fr	9.0 mm	170 mm	10

Endotracheal Tube

MURPHY EYE

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE



- Oral and nasal use
- Clear, flexible PVC
- Smooth Murphy Eye for alternative ventilation port
- Available cuffed, uncuffed, reinforced and with stylette
- Economical
- High volume, low pressure cuff
- Large barrel cuff
- Standard 15 mm ISO connector included with each tube



CUFFED

ITEM	ID	OD	LENGTH	SIZE	PK
1-7333-25	2.5 mm	3.5 mm	150 mm	10 Fr	10
1-7333-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7333-35	3.5 mm	4.9 mm	188 mm	14 Fr	10
1-7333-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7333-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7333-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7333-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7333-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7333-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7333-70	7.0 mm	9.6 mm	310 mm	28 Fr	10
1-7333-75	7.5 mm	10.2 mm	315 mm	30 Fr	10
1-7333-80	8.0 mm	10.9 mm	330 mm	32 Fr	10
1-7333-85	8.5 mm	11.5 mm	330 mm	34 Fr	10
1-7333-90	9.0 mm	12.1 mm	330 mm	36 Fr	10
1-7333-95	9.5 mm	12.7 mm	330 mm	38 Fr	10
1-7333-10	10.0 mm	13.6 mm	330 mm	40 Fr	10



CUFFED WITH STYLETTE

ITEM	ID	OD	LENGTH	SIZE	PK
1-7343-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7343-35	3.5 mm	4.9 mm	188 mm	14 Fr	10
1-7343-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7343-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7343-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7343-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7343-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7343-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7343-70	7.0 mm	9.6 mm	310 mm	28 Fr	10
1-7343-75	7.5 mm	10.2 mm	315 mm	30 Fr	10
1-7343-80	8.0 mm	10.9 mm	330 mm	32 Fr	10
1-7343-85	8.5 mm	11.5 mm	330 mm	34 Fr	10
1-7343-90	9.0 mm	12.1 mm	330 mm	36 Fr	10
1-7343-95	9.5 mm	12.7 mm	330 mm	38 Fr	10
1-7343-10	10.0 mm	13.6 mm	330 mm	40 Fr	10

UNCUFFED

ITEM	ID	OD	LENGTH	SIZE	PK
1-7330-20	2.0 mm	2.9 mm	145 mm	8 Fr	10
1-7330-25	2.5 mm	3.5 mm	150 mm	10 Fr	10
1-7330-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7330-35	3.5 mm	4.9 mm	185 mm	14 Fr	10
1-7330-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7330-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7330-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7330-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7330-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7330-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7330-70	7.0 mm	9.6 mm	310 mm	28 Fr	10

Tracheostomy Tube

TRACHEOSTOMY

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Obturator and tube clearly marked
- 15 mm Connector swivel adapter
- Soft PVC
- Secure neck band retainer
- Cushion neck band comfortable for patient
- Radiopaque strip embedded for exact location of tube
- Low pressure cuff

UNCUFFED

ITEM	ID	LENGTH	SIZE	PK
1-7390-30	3.0 mm	56.5 mm	12 Fr	10
1-7390-35	3.5 mm	56.5 mm	14 Fr	10
1-7390-40	4.0 mm	56.5 mm	16 Fr	10
1-7390-45	4.5 mm	56.5 mm	18 Fr	10
1-7390-50	5.0 mm	62 mm	20 Fr	10

CUFFED

ITEM	ID	LENGTH	SIZE	PK
1-7393-50	5.0 mm	62 mm	20 Fr	10
1-7393-55	5.5 mm	67 mm	22 Fr	10
1-7393-60	6.0 mm	67 mm	24 Fr	10
1-7393-65	6.5 mm	77 mm	26 Fr	10
1-7393-70	7.0 mm	80 mm	28 Fr	10
1-7393-75	7.5 mm	84 mm	30 Fr	10
1-7393-80	8.0 mm	90 mm	32 Fr	10
1-7393-85	8.5 mm	92 mm	34 Fr	10
1-7393-90	9.0 mm	97 mm	36 Fr	10
1-7393-95	9.5 mm	97 mm	38 Fr	10
1-7393-10	10.0 mm	97 mm	40 Fr	10

Introducer

BOUGIE-TO-GO

STERILE | LATEX FREE | SINGLE USE | DISPOSABLE

SunMed



- Innovative compact packaging designed specifically for use by EMS professionals
- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Fits in 6 mm to 11 mm tubes


ITEM	DESCRIPTION	SIZE	PK
9-0212-82	Adult, Coudé Tip	15 Fr x 60 cm	10



The Bougie-To-Go is specially packaged to fit in smaller compartments, bags, and even pockets, making it ideal for first responders

Intubation
INTUBATION KIT

LATEX FREE



Complete intubation kits available with conventional, fiber optic or fiber optic disposable laryngoscope blades

Each kit includes:

- 8 Oropharyngeal color coded Guedel airways (sizes 40 mm-100 mm)
 - 1 Penlite chrome plated laryngoscope handle
 - 1 Medium chrome plated laryngoscope handle
 - 5 Miller laryngoscope blades (sizes 0-4)
 - 4 Macintosh laryngoscope blades (sizes 1-4)
 - 4 Extra laryngoscope lamps (2 small, 2 large)
 - 1 Adult Magill forceps
 - 1 Pediatric Magill forceps
 - 2 "C" batteries
 - 2 "AA" batteries
 - 4 Packets of sterile lubricating jelly
 - 2 Endotracheal uncuffed tubes (sizes 2.5, 4 mm)
 - 8 Endotracheal cuffed tubes (sizes 5-10 mm)
- 4 Nasopharyngeal airways (sizes 20, 24, 28, 32 Fr)
 - 1 10 mL syringe
 - 1 Endotracheal tube introducer (bougie) with coude tip
 - 1 Adult stylette
 - 1 Nylon airway management case

ITEM	DESCRIPTION	PK
8-1050-95	Fiber Optic Intubation Kit Includes items above + Resuable FO Laryngoscope Blades	1
8-1050-85	Conventional Intubation Kit Includes items above + Conventional Laryngoscope Blades	1
8-1050-94	Single-Use Fiber Optic Intubation Kit Includes items above + Disposable FO Single Use Laryngoscope Blades	1
8-1050-80	Nylon Case Only	1

GreenLine® Fiber Optic Laryngoscope

DISPOSABLE LARYNGOSCOPE KITS


LATEX FREE | SINGLE USE | DISPOSABLE

SunMed




- Ideal for ambulances, paramedics and first responders
- Variety of blade/handle combinations in most popular sizes
- All blades are single patient use and disposable
- Available with either single use or reusable handles
- Components arrive in original packaging inside sturdy canvas pouch

ITEM	KIT INCLUDES	PK
5-5332-55	2 Disposable Standard Handles (5-5338-12) Macintosh Blades: 2 (5-5332-02), 3 (5-5332-03), 4 (5-5332-04) Miller Blades: 2 (5-5333-02), 3 (5-5333-03) Canvas Pouch	1
5-5332-57	Disposable Standard Handle (5-5338-12) Disposable Stubby Handle (5-5388-11) Macintosh Blades: 2 (5-5332-02), 3 (5-5332-03), 4 (5-5332-04) Miller Blades: 2 (5-5333-02), 3 (5-5333-03) Canvas Pouch	1
5-5336-54	2 Disposable Penlite Handles (5-5338-10) All-Metal Macintosh Blade: 1 (5-5335-01) All-Metal Miller Blades: 00 (5-5336-20), 0 (5-5336-00), 1 (5-5336-01) Canvas Pouch	1
5-5333-58	2 Disposable Standard Handles (5-5338-12) Macintosh Blades: 2 (5-5332-02), 3 (5-5332-03), 4 (5-5332-04), Miller Blades: 0 (5-5333-20), 1 (5-5333-01), 2 (5-5333-02), 3 (5-5333-03), Canvas Pouch	1
5-5333-57	Chrome Plated Med. Handle (5-0236-09) Macintosh Blades: 2 (5-5332-02), 3 (5-5332-03), 4 (5-5332-04) Miller Blades: 0 (5-5333-20), 1 (5-5333-01), 2 (5-5333-02), 3 (5-5333-03) Canvas Pouch	1
5-5335-47	Chrome Plated Med. Handle (5-0236-09) All-Metal Macintosh Blade: 2 (5-5335-02), 3 (5-5335-03), 4 (5-5335-04) All-Metal Miller Blades: 0 (5-5336-00), 1 (5-5336-01), 2 (5-5336-02), 3 (5336-03) Canvas Pouch	1




5-5332-55




5-5332-57




5-5336-54



5-5336-58



5-5333-57



5-5335-47


Individual blades and handles are also available for purchase.

Manual Resuscitator / Ventilator

AIRFLOW™

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ventlab



AirFlow™ MANUAL RESUSCITATOR

- Single patient use resuscitation/ventilation bag available in five sizes/volumes
- All series include bacterial/viral filter, oxygen reservoir and mask
- Available with or without integrated pressure manometer
- Integrated color-coded pressure manometer effectively promotes proper pressure delivery (<20 cm H₂O); the tri-color fields alert to the risk level of delivered pressure
- Accurate pressure delivery helps reduce gastric insufflation, thus preventing vomiting, aspiration and resulting pneumonia
- Face mask proportionately sized to resuscitation bag size
- Select from a variety of components to assemble a custom resuscitation bag to meet your unique patient needs

AIRFLOW™ + INTEGRATED PRESSURE MANOMETER

ITEM	SIZE	VOLUME	STROKE VOLUME	BODY MASS	PK
AF3140MB	Infant	300 mL	150-170 mL	≤ 10 kg (22 lbs)	10
AF4140MB	Child	550 mL	330-350 mL	≤ 22 kg (48 lbs)	10
AF2140MB	Pediatric	1000 mL	600-620 mL	≤ 40 kg (88 lbs)	10
AF5140MB	Adult	1500 mL	640-660 mL	> 40 kg (88 lbs)	10
AF1140MB	Adult	1900 mL	950-960 mL	> 40 kg (88 lbs)	10

AIRFLOW™ STANDARD (WITHOUT PRESSURE MANOMETER)

ITEM	SIZE	VOLUME	STROKE VOLUME	BODY MASS	PK
AF3040MB	Infant	300 mL	150-170 mL	≤ 10 kg (22 lbs)	10
AF4040MB	Child	550 mL	330-350 mL	≤ 22 kg (48 lbs)	10
AF2040MB	Pediatric	1000 mL	600-620 mL	≤ 40 kg (88 lbs)	10
AF5040MB	Adult	1500 mL	640-660 mL	> 40 kg (88 lbs)	10
AF1040MB	Adult	1900 mL	950-960 mL	> 40 kg (88 lbs)	10

CUSTOMIZE YOUR RESUSCITATION / VENTILATION BAG

Choose from the following components to create a custom resuscitation/ventilation bag to best meet your needs:

OXYGEN RESERVOIR



INFLATABLE O₂ BAG

Closed system oxygen reservoir delivers higher FDO₂ and provides visual indication of proper airflow



CORRUGATED TUBING

For use with secondary oxygen supply; available in 24" and 48" lengths



COLLAPSIBLE/EXTENDABLE O₂ ACCUMULATOR

Available only for AF1000 series

CONVENIENCE COMPONENTS



FACE MASKS

Inflatable and non-inflatable face masks available



FLEX NECK

Removable, flexible joint allows corrugated tubing to be joined with the patient valve; also reduces the torque on ET tube



HAND STRAP

Maintain a positive grip with a hand strap
*Only available on inflatable bag reservoir



SAFETY COMPONENTS



INTEGRATED PRESSURE MANOMETER

Ventlab's innovative color-coded dial clearly denotes pressure delivery:

Green - Target pressure level <20cm H₂O
Yellow - Potentially dangerous pressure level
Red - Dangerous pressure level



POP-OFF VALVES

40 cm/H₂O & 25 cm/H₂O
Allows clinician to use an appropriate pop-off level



FILTERS

A bacterial/viral filter provides effective protection for the healthcare provider during expiration



CO₂ INDICATOR

Choose from: STAT-Check II, MaxCap, or FLOCAP for effective assessment of CO₂ exchange



PEEP VALVES


Adjustable PEEP valves;
5-20 cm H₂O

Manual Resuscitator / Ventilator

SAFET™ T-PIECE RESUSCITATOR

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

SunMed



Integrated, color-coded pressure manometer for effective monitoring of ventilatory pressures

Variety of mask sizes available to meet patient and hospital needs

PEEP dial allows clinician to set PEEP; ventilation is provided by occluding PEEP dial hole intermittently at the desired respiratory rate allowing for inspiration and exhalation to occur

Expandable circuit that allows up to 18" in length


Standard pressure relief valve designed to lower the risk of over-pressurization

Oxygen tubing with red Fits-All universal connector

Inspiratory pressure controller allows clinician to set the PIP

Side-by-Side Comparison

Current NRP guidelines recommend the use of either self-inflating, flow-inflating or T-piece resuscitators for manual ventilation of neonates. This side-by-side comparison highlights the additional features of a T-piece resuscitator compared to a BVM and hyperinflation system.




Neonatal Resuscitator Features + Functionality	SafeT™	BVM	Hyperinflation
Continuous O ₂ Flow	●		●
Pre-Set + Deliver Consistent PIP	●		
Pre-Set + Deliver Consistent PEEP	●	●	●
Can Provide CPAP	●		●
Pre-set Maximum Pressure + Adjustability During Resus.	●		
Integrated Pressure Manometer	●	●	●
Pressure Relief Valve	●	●	●
Hand Compression Required		●	●
Portable	●	●	●



- Allows clinician to deliver consistent PIP and PEEP pressures
- PIP and PEEP ensures greater safety when ventilating babies' fragile lungs, particularly in an emergency situation
- Integrated pressure manometer provides clear visualization of ventilatory pressure ranges
- Pressure relief valve prevents over-pressurization
- Includes a cushioned face mask for proper sealing during mask ventilation
- Clear controller and manometer dials are easy to see, even in dim lighting
- Universal red Fits-All O₂ tubing connector allows immediate visualization of proper oxygen connection

ITEM	DESCRIPTION	MASK	CA
ST700	SafeT™ T-Piece Resuscitator	X-Small Neonate	10
ST701	SafeT™ T-Piece Resuscitator	Small Neonate	10
ST702	SafeT™ T-Piece Resuscitator	Small Infant	10
ST703	SafeT™ T-Piece Resuscitator	Infant	10
ST704	SafeT™ T-Piece Resuscitator	Neonate Round	10
ST705	SafeT™ T-Piece Resuscitator	Infant Round	10
ST707	SafeT™ T-Piece Resuscitator	NO MASK	10




Manual Resuscitator / Ventilator

SAFE SPOT™

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ventlab



SafeSpot™

- Single patient use resuscitation/ventilation system designed specifically for infant patients
- Bag features ultra-thin ‘spots’ to provide an enhanced tactile feel for lung compliance and increased patient safety
- Integrated pressure manometer promotes proper delivery (<20 cm H₂O) to guard against over-pressurization of babies’ fragile lungs; the tri-color fields alert to the risk level of delivered pressure
- Accurate pressure delivery helps reduce barotrauma or pneumothoraces
- Pop-off valve (40 cm H₂O and 25 cm H₂O available) enables clinicians to use the pop off level appropriate for a given patient population creating a greater margin of safety for the baby
- Optional Vent Mask in various sizes available
- Select from a variety of components to assemble a custom resuscitation bag to meet your unique patient needs

ITEM	SIZE	VOLUME	STROKE VOLUME	BODY MASS	PK
SS3200MB	Infant	250 mL	150-170 mL	≤ 10 kg (22 lbs)	10

CUSTOMIZE YOUR MANUAL RESUSCITATOR

Choose From the following components to create a custom resuscitation/ventilation bag to best meet your needs:

SAFETY COMPONENTS



INTEGRATED PRESSURE MANOMETER

Ventlab's innovative color-coded dial clearly denotes pressure delivery:

- Green** - Target pressure level <20 cm H₂O
- Yellow** - Potentially dangerous pressure level
- Red** - Dangerous pressure level



POP-OFF VALVES

40 cm H₂O & 25 cm H₂O
Allows clinician to use an appropriate pop-off level



FILTERS

A bacterial/viral filter provides effective protection for the healthcare provider during expiration



CO₂ INDICATOR

Choose From STAT-Check II or MaxCap for effective assessment of CO₂ exchange



PEEP VALVES

Adjustable PEEP valves;
5-20 cm/H₂O



SafeSpot™, SS3200-Series
Shown with integrated pressure manometer, pop-off valve, inflatable mask and corrugated tubing

SafeSpot™ bags come standard with an inflatable Vent Mask

OXYGEN RESERVOIRS



INFLATABLE O₂ BAG

Closed system oxygen reservoir delivers higher FDO₂ and provides visual indication of proper airflow



CORRUGATED TUBING

For use with secondary oxygen supply; available in 24" and 48" lengths

CONVENIENCE COMPONENTS



FACE MASKS

Inflatable and non-inflatable face masks available



FLEX NECK

Removable, flexible joint allows corrugated tubing to be joined with the patient valve; also reduces the torque on ET tube

Gastric Lavage

TUM-E-VAC® GASTRIC LAVAGE KIT

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox

Tum-E-Vac® Gastric Lavage Kits are fully assembled, disposable kits that promote quick, efficient lavage of stomach contents. This completely closed system is designed to protect healthcare professionals from potential exposure to stomach contents and charcoal solutions during gastric lavage or decontamination procedures.



- Pre-assembled kits are quick and easy to use, reducing the time for lavage procedures
- Promotes quick, efficient lavage of stomach contents
- Closed system helps prevent exposure to stomach contents and charcoal solutions
- Easy-to-fill top port
- Available with Universal Solution Spike (for use with pre-filled solution containers)
- No clean up or disinfecting procedures necessary
- All kits come standard with Bite Block

ITEM	DESCRIPTION	PK
2054	Tum-E-Vac® + 3 L Solution Bag + 32 Fr Lavage Tubing	10
2055	Tum-E-Vac® + 3 L Solution Bag	10
2075	Tum-E-Vac® + Universal Solution Spike	10
2060	Tum-E-Vac® Bite Block Only	50


Gastric Lavage

LAVAGE TUBING

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox

Radiopaque strip
confirms
proper tube
placement



Reference marks at
18", 20" and 22" From
distal end



- Open bullet tip design with multiple eyelets to create free-flowing movement with low resistance
- Soft clear tube with smooth tip
- Suction connector on proximal end
- Radiopaque stripe
- All tubing is 48" long with reference marks at 18", 20" and 22" from distal end
- Individually packaged

ITEM	DESCRIPTION	SIZE	LENGTH	PK
2124	Lavage Tubing	24 Fr	48"	10
2128	Lavage Tubing	28 Fr	48"	10
2132	Lavage Tubing	32 Fr	48"	10
2136	Lavage Tubing	36 Fr	48"	10

Airway Accessory

JAW SPREADER KIT

LATEX FREE | SINGLE USE | DISPOSABLE

SunMed

- Helps facilitate breathing during convulsions or seizures
- Universal design firmly holds the tongue down while keeping jaw open in both children and adults
- Dependable, durable construction – will not break or splinter even under extreme pressure

ITEM	DESCRIPTION	PK
9-3042-00	Kit contains: Jaw Spreader (fits children and adults); Pediatric Berman Airway (6 cm); Adult Berman Airway (9 cm); Case	1

Patient Safety

BITE BLOCK

LATEX FREE | SINGLE USE | DISPOSABLE

Ethox

- Assists with naso-gastric tube placement
- Accommodates up to a 40 Fr tube
- Individually packaged


ITEM	DESCRIPTION	PK
2060	Bite Block	50

Wound Irrigation


COMBIPORT® + COMBIGUARD II™

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

Ethox



Working together, a saline bag uses a Combiport® attached to a Combiguard II™ to effectively cleanse the wound without splash-back on the caregiver



COMBIPORT® WOUND IRRIGATION DEVICE

- Attaches to a standard saline bag to administer a high volume of fluid in a short period of time
- Irrigation helps remove foreign bodies, reduce pathogen numbers and moisten wound tissue
- Patented spike attaches to standard 500 mL or 1000 mL saline bag
- For use only with Combiguard Irrigation Splash Guard (not sold separately)
- Removable cap for easy stop/start during cleansing
- Easy to assemble and use

STERILE

ITEM	DESCRIPTION	PK
10222	Combiguard II™ Irrigation Splash Guard	50
10100	Combiguard II™ + Combiport®	50

NON-STERILE

ITEM	DESCRIPTION	PK
10222NS	Combiguard II™ Irrigation Splash Guard Non-Sterile	500

COMBIGUARD II™ IRRIGATION SPLASH GUARD

- Egg-shaped cup covers wound site to prevent irrigation splash-back, helping to protect clinician
- Dual outlet ports enable contaminated fluids to safely escape
- Fluid delivery inlet delivers a stream equivalent to an 18 g needle
- Angled delivery inlet directs fluid away from cleansed area toward outlet ports
- Rounded edges provide a smooth surface against the skin
- Fits on 30-40 cc luer-tipped syringe
- Can be combined with Combiport Wound Irrigation Device


Rapid Infusion

INFU-SURG® PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Ethox**INFU-SURG®**

- Ideal for intravenous and intra-arterial infusion treatment including A-line pressure monitoring
- Pressure relief valve prevents over inflation (330 mm Hg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Reinforced hanging loop allows for quick and easy hanging on IV pole
- Color-coded gauge makes for accurate pressure monitoring (150-300 mm Hg)
- Three-way stopcock ensures precise control of pressure
- Piston-style gauge allows for 360° visibility



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation


ITEM	DESCRIPTION	VOLUME	PK
4005H	Infu-Surg® with Hook	500 mL	5
4010H	Infu-Surg® with Hook	1000 mL	5
4030	Infu-Surg®	3000 mL	5
4030H	Infu-Surg® with Hanger Loop	3000 mL	5

Forceps

OPEN MAGILL FORCEPS

LATEX FREE

SunMed



- Surgical stainless steel

ITEM	DESCRIPTION	SIZE	LENGTH	PK
9-3016-02	Magill Forceps	Child	4.5"	1
9-3016-03	Magill Forceps	Adult	5.5"	1

Instruments

TONGUE SEIZING FORCEPS

LATEX

SunMed



- Surgical stainless steel
- Includes latex rubber tips

ITEM	DESCRIPTION	SIZE	PK
9-3021-50	Tongue Seizing Forceps	3"	1
9-3021-52	Latex Replacement Tips	--	1 pair

Forceps

CLOSED MAGILL FORCEPS

LATEX FREE

SunMed



- Surgical stainless steel

ITEM	DESCRIPTION	SIZE	LENGTH	PK
9-3018-01	Closed Tip Magill Forceps	Infant	3"	1
9-3018-02	Closed Tip Magill Forceps	Child	4"	1
9-3018-03	Closed Tip Magill Forceps	Adult	5"	1

Mask

CPR POCKET MASK

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

SunMed



- Designed with breathing filter and one-way valve

ITEM	DESCRIPTION	PK
8-5070-03	CPR Pocket Mask with Large Filter and One-way Valve	1

Stethoscope SPRAGUE

LATEX FREE

SunMed



- Three different size bells eliminates the need for multiple stethoscopes for infants, children and adults
- Various heads for acoustical responses in low, medium and high pitch sounds
- Heavy wall 22" tubing and choice of ear tips help reduce sound leakage and extraneous noise
- Complete stethoscope kit includes 1 pair black plastic ear tips, 3 open bells (infant, child, adult), 1 pair soft rubber (mushroom type) ear tips, adult and pediatric diaphragms
- 29" Overall length

ITEM	DESCRIPTION	PK
6-0028-00	Sprague Rappaport Stethoscope	1

Stethoscope DUAL HEAD

LATEX FREE

SunMed



- Molded one-piece Y tubing
- High sensitivity - lightweight diaphragm
- Non-chill ring
- Anodized aluminum chestpiece
- Chrome plated brass binaurals
- White plastic eartips
- Superior acoustics
- 31" Overall length

ITEM	DESCRIPTION	COLOR	PK
6-0027-40	Dual Head Stethoscope	Black	1
6-0027-41	Dual Head Stethoscope	Red	1
6-0027-42	Dual Head Stethoscope	Grey	1
6-0027-43	Dual Head Stethoscope	Green	1