

Flexi IV-Boards™

Our Flexi IV-Boards™ have been specifically designed for patient comfort and stability of the limb when securing IV lines. The firm yet flexible aluminium insert allows the Flexi IV-Boards™ to be moulded and hold its shape on the patient.

FEATURES AND BENEFITS

- Soft and pliable
- Reduce discomfort to limb
- Easy to apply and remove
- Available in 7 sizes (Micro Neonatal to Adult)
- Firm, yet flexible aluminium insert to help keep profile
- Latex free

REF	Description
APM601	Flexi IV-Boards™ Micro Neonatal Armboard.
APM602	Flexi IV-Boards™ Neonatal Armboard.
APM603	Flexi IV-Boards™ Infant Armboard.
APM603.5	Flexi IV-Boards™ Small Paediatric Armboard.
APM604	Flexi IV-Boards™ Medium Paediatric Armboard.
APM604.5	Flexi IV-Boards™ Large Paediatric Armboard.
APM605	Flexi IV-Boards™ Adult Armboard.

Micro-Neonatal

APM601

Neonatal

APM602

Infant

APM603

Small Paediatric

APM603.5

Medium Paediatric

APM604

Large Paediatric

APM604.5

Adult

APM605

* Available with securing straps

* Product subject to change without notice



Mediplas Respiratory Products Sdn Bhd
No. 7 Jalan KIP 3, Taman Perindustrian KIP 52200 Kuala Lumpur, Malaysia
Phone: +60 3 6272 7625 Fax: +60 3 6273 6452
E-mail: sales@mediplas.com.my

EC REP Obelis s.a
Boulevard Général Wahis 53, 1030 Brussels, BELGIUM
Phone: +(32) 2 732 5954 Fax: +(32) 2 732 6003 E-Mail: mail@obelis.net

